

## CLL + CNV Panel

### Features

#### **Unparalleled uniformity and high depth of coverage**

- Detect low-frequency SNVs and indels with confidence

#### **CNV detection ranging from loss of single exons to full chromosome arms and trisomy 12**

- Profile your samples for CNVs in the 5 most commonly aberrant regions in CLL

#### **Time savings**

- Replace multiple assays with a single NGS panel, increasing throughput and reducing turnaround time

#### **Complimentary data analysis software**

- Analyse your data with Interpret, OGT's powerful and easy-to-use analysis solution for accurate identification of all variants and CNVs

### Introduction

Chronic lymphocytic leukaemia (CLL) is the most common type of leukaemia in adults. A wide variety of chromosomal abnormalities are associated with CLL, ranging from single nucleotide variants (SNVs) and insertions/deletions (indels) up to large copy number variations (CNVs), including trisomies.

The SureSeq™ CLL + CNV Panel has been designed in collaboration with recognised cancer experts to detect 13 key genes and 5 chromosomal regions implicated in CLL progression (Table 1). The SureSeq CLL + CNV Panel alleviates the burden of running multiple assays and streamlines your CLL research to deliver a comprehensive genomic profile for each CLL sample using a single workflow.

### Contains the latest evidence-based content

Investigating both chromosomal aberrations and SNVs/indels is imperative to advance research into CLL progression and treatment. Structural abnormalities are common in CLL and found in more than 80% of CLL cases, the most frequent being del(13q), del(11q), del(17p), del(6q) and trisomy 12<sup>1</sup>. Some of these CNVs cover important tumour suppressors, such as del(17p) resulting in the loss of the *TP53* gene. More recently, other genes have also been found to be mutated in CLL, including *NOTCH1*, *SF3B1*, *MYD88* and *BIRC3*, adding to the genomic complexity of this leukaemia<sup>2</sup>.

Due to this genetic heterogeneity, current analysis strategies for CLL require multiple methods to obtain a comprehensive genetic picture, often using microarray or fluorescence *in situ* hybridisation (FISH) to detect structural abnormalities in combination with NGS for somatic variants. With OGT's SureSeq CLL + CNV Panel, you can now obtain a more complete understanding of the genetic makeup of CLL progression in each sample using a single assay.

Genes	<i>ATM</i> , <i>BIRC3</i> , <i>BRAF</i> , <i>BTK</i> , <i>CXCR4</i> , <i>KRAS</i> , <i>MYD88</i> , <i>NOTCH1</i> , <i>PLCG2</i> , <i>SAMHD1</i> , <i>SF3B1</i> , <i>TP53</i> and <i>XPOT</i>
CNVs	17p (covering <i>TP53</i> ), 11q (covering <i>ATM</i> ), 13q (covering <i>RB1/DLEU2/DLEU7</i> ), 6q (6q23.2-6q23.3 covering <i>MYB</i> ) and Trisomy 12.
Sample tracking	<i>SRY</i> + 24 SNP profiling panel <sup>3</sup>

Table 1: The SureSeq CLL + CNV Panel targets the 5 most common chromosomal regions implicated in CLL and 14 genes, including the *SRY* gene and 24 SNPs for easy sample tracking.

### Superior coverage uniformity allowing reliable variant and somatic CNV detection

OGT's expert bait design delivers outstanding uniformity and depth of coverage, capable of detecting low frequency SNVs and indels down to 1% low frequency SNVs and indels down to 1% variant allele frequency (VAF) in 14 genes (Figure 1), including the *SRY* gene and 24 SNPs to allow for easy sample tracking<sup>3</sup>.





Figure 1: Illustration of the excellent uniformity and high depth of coverage allowing confident detection of **A** a *SF3B1* exon 15 hotspot variant Lys700Glu with 4.8% allele frequency and **B** a *TP53* exon 4 frameshift deletion (*TP53* c.124del) with frequency 38.9%.

The SureSeq CLL + CNV Panel covers the 5 most common CNVs in CLL with the capability of detecting VAF down to 10% , corresponding to 20% tumour content. Compared to array data, often considered the gold standard for CNV detection, the events reported with the SureSeq CLL + CNV Panel were 100% concordant, even in genomic regions containing multiple aberrations (Figures 2 – 3). More so, facilitated by OGT’s excellent bait design, loss-of-heterozygosity (LOH) can be identified. With a CNV size detection range from single exon to whole gene, up to complete loss of a chromosomal arm and whole chromosome gains (trisomy 12), your data provides a more comprehensive genetic picture for each sample from a single assay.

### NGS



### Array





Figure 2: 42.7 Mb deletion of 11q covering *ATM*.

### NGS



### Array




Figure 3: 0.6 Mb biallelic loss called within a larger ~1 Mb single allele deletion in the region covering *DLEU2/DLEU1/DLEU7* on chromosome 13q.

## Complimentary Interpret software

Interpret is OGT's powerful and easy-to-use data analysis solution, facilitating analysis and visualisation of a wide range of somatic variants and structural aberrations. Designed to work seamlessly with all SureSeq panels, Interpret perfectly complements the SureSeq CLL + CNV Panel, delivering fast and accurate detection of all SNVs, indels, LOH and CNVs covered by the panel. Following detection, all events can be readily visualised in the user-friendly variant browser, for an effortless translation of all your CLL data into meaningful results (Figure 4).




Figure 4: Following analysis, all variants and CNVs are visualised for easy interpretation in OGT's Interpret. In this example a trisomy 12 is detected, showing a reliable gain call across the whole chromosome.

### The OGT Partnership

Behind every sample is a life that can be improved through the right care decisions. The OGT partnership approach is key to providing the highest level of service, working closely with you to understand your unique challenges, customising our approach to meet your exact needs.

### Bespoke panel content

Does the SureSeq CLL + CNV Panel not meet your exact requirements? With OGT, you never have to sequence genes you're not interested in and can always modify each panel to what's most relevant for your research. Choose from our regularly updated, expert-curated library of pre-optimised cancer content to create your ideal custom SureSeq myPanel™ CLL Panel, or order the SureSeq CLL + CNV Panel right off the shelf.

### The SureSeq CLL + CNV Panel in numbers

Feature	Specification
Number of genes	14
Panel size	117 kb
Uniformity Coverage	>95% of the panel achieved >20% of the mean target coverage
Mean target coverage	>1000x
DNA input recommended	>500ng high quality DNA
Gene list	<i>ATM, BIRC3, BRAF, BTK, CXCR4, KRAS, MYD88, NOTCH1, PLCG2, SAMHD1, SF3B1, TP53</i> and <i>XPO1</i>
CNV list	17p (covering <i>TP53</i> ), 11q (covering <i>ATM</i> ), 13q (covering <i>RB1/DLEU2/DLEU7</i> ), 6q (6q23.2-6q23.3 covering <i>MYB</i> ) and Trisomy 12
Sample tracking	<i>SRY</i> + 24 SNP profiling panel <sup>3</sup>
Limit of Detection	SNVs/indels: capable of VAF of 1% within the 14 genes CNVs: capable of detecting VAF of 10% within the 5 chromosomal regions
CNV detection size	11q: - single exon to whole gene of <i>ATM</i> - > 5-10 Mb for the rest of the 11q arm
	17p: - single exon to whole gene of <i>TP53</i> - > 5-10 Mb for the rest of the 17p arm
	13q: - del(13)(q14) type I (short) and del(13)(q14) type II (larger) events covering <i>RB1/DLEU2/DLEU7</i> - > 10-20 Mb for the rest of the 13q arm
	6q (6q23.2-6q23.3): - single exon to whole gene of <i>MYB</i> plus 1 Mb flanking sequence on either side
	Trisomy 12: - whole chromosome
LOH detection size	11q and 17p: 5-10 Mb 13q: 10-20 Mb
Samples per MiSeq® v2 run	16 samples/ run

### Ordering information

UK +44 (0) 1865 856800

US +1 914 467 5285

contact@ogt.com

ogt.com

Product	Contents	Cat. No.
SureSeq CLL + CNV Complete NGS Workflow Solution (24 reactions)	Bundle of Enrichment baits sufficient for 3 pools of 8 samples; 1x Universal NGS Library Preparation Kit (24), containing PCR primers and enzymes, 1x Universal NGS Index Adapter Kit (24), 1x Universal NGS Hyb & Wash Kit (24), 1x Universal NGS Bead Kit (24)	780004-24
SureSeq CLL + CNV Complete NGS Workflow Solution (96 reactions)	Bundle of Enrichment baits sufficient for 12 pools of 8 samples; 1x Universal NGS Library Preparation Kit (96), containing PCR primers and enzymes, 1x Universal NGS Index Adapter Kit (96), 1x Universal NGS Hyb & Wash Kit (96), 1x Universal NGS Bead Kit (96)	780004-96
SureSeq CLL + CNV Panel (24 reactions)	Enrichment baits sufficient for 3 pools of 8 samples; Interpret Software	770004-24
SureSeq CLL + CNV Panel (96 reactions)	Enrichment baits sufficient for 12 pools of 8 samples; Interpret Software	770004-96
Universal NGS Workflow Solution (24)	Bundle of 1x Universal NGS Library Preparation Kit (24), containing PCR primers and enzymes, 1x Universal NGS Index Adapter Kit (24), 1x Universal NGS Hyb & Wash Kit (24), 1x Universal NGS Bead Kit (24)	770500-24
Universal NGS Workflow Solution (96)	Bundle of 1x Universal NGS Library Preparation Kit (96), containing PCR primers and enzymes, 1x Universal NGS Index Adapter Kit (96), 1x Universal NGS Hyb & Wash Kit (96), 1x Universal NGS Bead Kit (96)	770500-96

### References

1. Döhner *et al.*, *N Engl J Med* 2000;343:1910–1916
2. Rossi *et al.*, *Blood* 2013;121:1403–1412
3. Pengelly *et al.*, *Genome Med* 2013;5:89



A Sysmex Group Company

**What binds us,  
makes us.**

Oxford Gene Technology Ltd., Begbroke Science Park, Woodstock Road, Begbroke, Oxfordshire, OX5 1PF, UK

**SureSeq:** For Research Use Only; Not for Diagnostic Procedures. This document and its contents are © Oxford Gene Technology IP Limited – 2022. All rights reserved. OGT™ and SureSeq™ are trademarks of Oxford Gene Technology IP Limited. Laboratories must undertake appropriate validation as outlined in local Directives and Guidelines.

