A Sysmex Group Company

## Instructions For Use

REF: DES500L, DFS500L, DSS500L, DES1000L

### DAPI Antifade

PROFESSIONAL USE ONLY

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Further information available at [www.cytocell.com](http://www.cytocell.com)

#### Intended purpose

DAPI Antifade solution is intended to aid visualisation of hybridised DNA FISH probes using fluorescent microscopy. For optimum results we recommend the use of DAPI Antifade ES with all Cytocell IVD FISH probes.

Table 1: DAPI Antifade Solution.

Catalogue Number	Description	Volume
DES500L	DAPI Antifade ES: 0.125µg/ml DAPI (4,6-diamidino-2-phénylindole)	500µl
DSS500L	DAPI Antifade SS: 0.0625µg/ml DAPI (4,6-diamidino-2-phénylindole)	500µl
DFS500L	DAPI Antifade FS: 1µg/ml DAPI (4,6-diamidino-2-phénylindole)	500µl
DES1000L	DAPI Antifade ES: 0.125µg/ml DAPI (4,6-diamidino-2-phénylindole)	1000µl

This DAPI Antifade kit contains only one of the products from the range of DAPI solutions noted above.

#### Materials Provided

DAPI Antifade solution: 500µl per vial or 1000µl per vial

#### Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- This product is intended to be used with Cytocell products or in an LDT (Laboratory Developed Test). If used as part of an LDT it is the responsibility of the laboratory developing the test to validate the test before use in a clinical setting.
- Wear gloves when handling Hybridisation solution, DNA probes and DAPI counterstain.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- Dispose of all hazardous materials according to your institution's guidelines for hazardous waste disposal.
- Operators must be capable of visually distinguishing between red, blue and green.
- Users should ensure that the correct concentration of DAPI Antifade is selected according to the products used.

#### Storage and Handling

Store the DAPI Antifade solution between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. Store the DAPI Antifade solution in the dark. The DAPI must be allowed to warm up at room temperature for 10 minutes prior to use.

#### Materials Necessary but not Supplied

##### Equipment

Variable volume micropipettes and tips range 1 µl - 200µl.

##### Fluorescence Microscope Recommendation

Microscope filters should be selected according to the products used. The DAPI filter should be used for observing DAPI Antifade. Ensure that a single filter is used for analysing blue probes. Avoid mixing DAPI Antifade with microscope immersion oil as this will obscure signals.

Check the fluorescence microscope before use to ensure it is operating correctly. Use immersion oil that is suitable for fluorescence microscopy and formulated for low autofluorescence. Follow manufacturers' recommendations in regards to the life of the lamp and the age of the filters.

#### Recommended FISH protocol using Cytocell liquid DNA FISH probes.

(Note: Ensure that exposure of DAPI to laboratory lights is limited at all times)

- After removing a slide from the post hybridisation washes remove the excess solution from the slide on the edge of a tissue.
- Ensure that the DAPI solution is uniformly mixed with a pipette.
- Apply DAPI Antifade to each end of the slide:
  - for a 24x24mm coverslip: use 10 - 20µl of DAPI Antifade
  - for a full slide: use 50µl of DAPI Antifade
- Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- View with a fluorescence microscope.

#### Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at 4°C.

#### Limitations

Failure to adhere to the protocol may affect the performance and results.

#### Additional Information

For additional product information please contact the Cytocell Technical Support Department.

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**Patents and Trademarks**
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
**Cytocell Ltd.**