Instructions For Use
REF: LPH 064-S/ LPH 064

FAST PML/RAR α Translocation, Dual Fusion Probe

PROFESSIONAL USE ONLY

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Further information available at www.ogt.com



Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Information

The fusion gene PML/RAR α (RARA) is created by the t(15;17)(q24.1;q21.2) translocation, found in 98% of AML M3 Acute hypergranular promyelocytic leukaemia and 9% of AML overall^{1,2,3}. Breakpoints in the PML gene vary between intron 3 and exon 7a, which is in contrast to the RARA breakpoint, which remains constant in intron 2. Variant translocations include partner breakpoints in 11q23, 5q32 and 11q13. The presence of three-way translocations indicate that the critical event occurs on the der(15), which always receives the distal end of chromosome 17. The PML protein is a transcription factor and RARA encodes a nuclear receptor. The fusion protein generated, PML-RARA, is a chimaeric transcription factor that operates as a dominant negative form of RARA and results in cells that are blocked at the promyelocytic stage of differentiation and then proliferate. Gain of function mutation allows repression of multiple genes and recruitment of DNA methyltransferases to promoters, allowing prolonged suppression. PML-RARA also activates components of the WNT signalling pathway, promoting stem cell renewal. Immediate treatment of PML/RARA positive patients is critical as intravascular coagulation causes early death in 10-40% of cases². Terminal differentiation is induced by the use of all-trans-retinoic-acid (ATRA), which reactivates the RARA gene and degrades the PML-RARA fusion product and 80-90% of treated patients achieve complete remission^{1,2,4}. Despite the efficacy of ATRA, the death rate remains around 15-20%, but 70% of patients will achieve 3 year survival². Additional abnormalities found with PML/RARA translocations include trisomy 8, seen in one third of cases, del(7q) and del(9q).

Probe Specification

PML, 15q24.1, Red
RAR α , 17q21.1-q21.2, Green

The PML probe mix, labelled in red, consists of a 151kb probe centromeric to the PML gene and a 174kb probe telomeric to the PML gene. The RAR α probe mix, labelled in green, consists of a 167kb probe centromeric to the RAR α (RARA) gene, including the CASC3 gene, and a 164kb probe, including the telomeric end of the RAR α gene as well as the TOP2A and IGFBP4 genes.

Materials Provided

Probe: 50 μ l per vial (5 tests) or 100 μ l per vial (10 tests)

Amount of red PML probe: 16-20ng/test

Amount of green RAR α probe: 48-60ng/test

The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

Counterstain: 150 μ l per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125 μ g/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe the fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The kit should be stored between -25 $^{\circ}$ C to -15 $^{\circ}$ C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80 $^{\circ}$ C).
2. Variable volume micropipettes and tips range 1 μ l - 200 μ l.
3. Water bath with accurate temperature control at 72 $^{\circ}$ C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37 $^{\circ}$ C incubator.
14. Rubber solution glue.

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all fluorophores and DAPI simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells or cultured bone marrow cells fixed in Carnoy's fixative, that should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

FAST FISH Protocol – 1 hour hybridisation

(Note: Please ensure that exposure of the probe to laboratory lights is limited at all times).

Slide preparation

1. Spot the cell sample onto a glass microscope slide. Allow to dry.
2. Immerse the slide in 2xSSC for 2 minutes at room temperature (RT) without agitation.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 minutes at RT.
4. Allow to dry.

Pre-Denaturation

6. Remove the probe from the freezer and allow it to warm to RT.
7. Ensure that the probe solution is uniformly mixed with a pipette.
8. Remove 10 μ l of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
9. Place the probe and the sample slide to prewarm on a 37 $^{\circ}$ C (+/- 1 $^{\circ}$ C) hotplate for 5 minutes.

- Spot 10µl of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

- Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2 minutes.

Hybridisation

- Place the slide in a humid, lightproof container at 37°C (+/- 1°C) for one hour.

Post-Hybridisation Washes

- Remove the coverslip and all traces of glue carefully.
- Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
- Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
- Drain the slide and apply 10µl of DAPI antifade onto each sample.
- Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- View with a fluorescence microscope.

Standard FISH Protocol – overnight hybridisation

(Note: Please ensure that exposure of the probe to laboratory lights is limited at all times).

Slide preparation

- Spot the cell sample onto a glass microscope slide. Allow to dry.
- Immerse the slide in 2xSSC for 2 minutes at room temperature (RT) without agitation.
- Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 minutes at RT.
- Allow to dry.

Pre-Denaturation

- Remove the probe from the freezer and allow it to warm to RT.
- Ensure that the probe solution is uniformly mixed with a pipette.
- Remove 10µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
- Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
- Spot 10µl of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

- Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2 minutes.

Hybridisation

- Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

- Remove the coverslip and all traces of glue carefully.
- Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
- Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
- Drain the slide and apply 10µl of DAPI antifade onto each sample.
- Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- View with a fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at/or below RT.

Procedural Recommendations

- Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
- Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by Cytocell Ltd.
- The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators as these temperatures are critical for optimum product performance.
- The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
- Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

Expected Results

In a normal cell these probes should appear as discrete red and green spots, one for each homologue (resulting in a 2R, 2G conformation). In a t(15;17)(q24.1;q21) patient there should be two yellow fusion signals in addition to the red and green signals of the normal chromosomes 15 and 17 respectively (1R, 1G, 2Y).

Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

Additional Information

For additional product information please contact the CytoCell Technical Support Department.

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Patents and Trademarks
CytoCell is a registered trademark of Cytocell Ltd.
3. Heim and Mitelman, Wiley-Liss, Inc. 1995