A Sysmex Group Company

Instructions For Use
REF: LPS 045-S / LPS 045

RET Breakpart Probe

PROFESSIONAL USE ONLY

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Further information available at www.ogt.com

Intended Use

The RET Breakpart Probe is intended to identify, via fluorescence in situ hybridisation (FISH), rearrangements of the RET region on chromosome 10q at location q11.21. It is intended to be used on formalin-fixed paraffin-embedded (FFPE) tissue sections or tissue microarrays (TMAs). This product is for professional use only and is intended to be an adjunct to Clinical Cytogenetics.

Background

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied to the assessment of solid tumour biopsies as well, which can provide important information for prediction of tumour progression. Current methodologies, namely immunohistochemistry or blotting can provide information at the gene expression level but, when carried out on tissue sections (either cryostat or paraffin embedded), FISH can provide information at a gene level, *in situ*, at the precise site within the tumour. This can reveal cell-to-cell heterogeneity and enable the detection of small clones of genetically distinct cells.


Probe Information

The RET (*ret proto-oncogene*) gene at 10q11 encodes for a transmembrane tyrosine kinase receptor involved in the control of cell differentiation, cell proliferation, and cell survival¹.

Rearrangements involving the RET gene are recognised recurrent abnormalities seen in 1-2% patients with lung adenocarcinomas, where it is seen fused with KIF5B^{2,3}, and papillary thyroid carcinoma where it is seen fused to a number of different partner genes including: CCDC6, PRKAR1A and NCOA4^{4,5}. The features of the proteins encoded by all types of RET fusion gene are similar to those of ALK: coiled-coil domains in the N-terminal fusion partners cause the RET domains to dimerize, resulting in activation of RET tyrosine kinase in the absence of ligands³.

Probe Specification

RET, 10q11.21, Green
RET, 10q11.21, Red

The RET Breakpart probe consists of two probes (92kb, 303kb), labelled in green, situated proximal to the RET gene and covering markers G17488 and D10S1100 and two probes (351kb, 112kb), labelled in red, situated distal to the RET gene and covering markers D10S1176 and D10S2221.

Materials Provided

Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of red RET probe: 68-118ng/test

Amount of green RET probe: 411-615ng/test

The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

Counterstain: 150µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. Dispose of all hazardous materials according to your institution's guidelines for hazardous waste disposal.
6. Operators must be capable of visually distinguishing between red, blue and green.

Storage and Handling

The kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench-top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.
15. Tissue Pretreatment Kit (LPS 100).

Fluorescence Microscope Recommendation

Use a 100-watt mercury lamp and plan apochromat objectives x63 or x100 for optimal visualisation. Use a triple bandpass filter DAPI/FITC/Texas Red for optimal visualisation of the green and red fluorophores and DAPI simultaneously. Alternatively for viewing red and green fluorophores use dual bandpass filter FITC/Texas Red.

Check the fluorescence microscope before use to ensure it is operating correctly. Use immersion oil that is suitable for fluorescence microscopy and formulated for low autofluorescence. Follow manufacturers' recommendations in regards to the life of the lamp and the age of the filters.

Sample Preparation

The kit is designed for use on Formalin Fixed Paraffin Embedded (FFPE) tissue sections or Tissue Microarrays (TMA) and should be prepared according to the laboratory or institution guidelines. Use 4µm - 6µm thick FFPE tissue sections for FISH.

Tissue Sample Pretreatment

Tissue sample pretreatment should be done according to the laboratory or institution guideline. Use the Tissue Pretreatment Kit (LPS 100) for optimal results.

FISH Protocol

(Note: Ensure that exposure of the probe to laboratory lights is limited at all times)

Pre-Denaturation

1. Remove the probe from the freezer and allow it to warm to room temperature (RT). Briefly centrifuge tubes before use
2. Ensure that the probe solution is uniformly mixed with a pipette.
3. Remove 10µl - 15µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
4. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
5. Spot 10µl - 15µl of probe mixture onto the sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

6. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 5 minutes.

Hybridisation

7. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

8. Remove the coverslip and all traces of glue carefully.
9. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
10. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
11. Drain the slide and apply 10µl - 15µl of DAPI antifade onto each sample.
12. Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
13. View with a fluorescence microscope.

Comments

Hybridisation efficiency and tissue morphology are usually negatively correlated. Aggressive pretreatment procedures improving hybridisation efficiency (e.g. an extended enzyme digestion time) tend to destroy cell structure and tissue

morphology. However, mild pretreatment saving tissue structures may not be sufficient for probe penetration and successful FISH results.

The optimal length of heat pretreatment and enzyme digestion time will depend on the age of the block, the tissue composition, and quality of tissue fixation. Decrease enzyme digestion for core biopsies and any sections that contain few tumour cells or have large areas of necrosis. Handle these samples particularly carefully to avoid over-digestion.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark below 4°C.

Procedural Recommendations

1. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCell Ltd.
2. Use a calibrated thermometer for measuring temperatures of solutions, waterbaths and incubators as these temperatures are critical for optimum product performance.
3. The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
4. Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.
5. Over hybridisation can result in additional or unexpected signals.
6. Users should optimise the protocol for their own samples prior to using the test for diagnostic purposes.

Expected Results

In a normal cell, two red/green fusion signals are expected (2F). In a cell with a balanced RET rearrangement, the expected signal pattern will be 1 fusion, 1 red and 1 green (1R, 1G, 1F). Other signal patterns are possible in aneuploid/unbalanced specimens.

Known Cross-Reactivity

No known cross-reactivity.

Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

Failure to adhere to the protocol may affect the performance and results.

This kit has not been validated for purposes outside of the intended use stated.

Additional Information

For additional product information please contact the CytoCell Technical Support Department.

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Patents and Trademarks
CytoCell is a registered trademark of CytoCell Ltd.