A Sysmex Group Company

Instructions For Use
REF: LPU 016-S / LPU 016

Kallmann (KAL1)/STS Probe Combination

PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com


Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Information

Kallmann syndrome (KS) is a developmental disease characterised by olfactory deficiency³ and hypogonadotropic hypogonadism (HH), which is responsible for the absence of spontaneous puberty¹. Kallmann syndrome is a heterogeneous developmental genetic disorder affecting approximately 1 in 8,000 males and 1 in 40,000 females². Reports indicate three modes of inheritance: X-linked, autosomal dominant and autosomal recessive^{1,4}. It has been shown that mutations in the KAL1 gene on Xp22.3 result in the X-linked form of KS⁵. KAL1 consists of 14 exons and extends over 200kb and abnormalities of KAL1 reported in patients with KS include missense and nonsense mutations, splice site mutations, intragenic deletions and submicroscopic chromosomal deletions involving the entire KAL1 gene⁷. Steroid Sulphatase Deficiency (STS) (also known as X-linked Ichthyosis)⁸ is the second most common type of ichthyosis and one of the most frequent human enzyme deficiency disorders. Deficiency of the STS enzyme is known to be responsible for dark, adhesive and regular scaling of the skin⁹. The gene encoding this protein maps to the distal short arm of chromosome X, which escapes X-chromosome inactivation and has the highest ratio of chromosomal deletions among all genetic disorders¹⁰. Complete deletions of the STS gene have been found in more than 90% of patients¹¹. The deletions can extend to involve neighbouring genes, causing contiguous gene defects. Therefore, STS may be associated with KS¹².

Probe Specification

KAL1, Xp22.31, Red
STS, Xp22.31, Green
DXZ1, Xp11.1-q11.1, Green

The KAL1 probe is 334kb, labelled in red and covers the entire KAL1 gene and the DXS278 and DXS7053 markers. The STS probe is 282kb, labelled in green and covers the STS gene and most of the HDHD1A and STS genes. The probe mix also contains a control probe for the X centromere (DXZ1), labelled in green.

Materials Provided

Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)
Amount of red KAL1 probe: 10-12ng/test
Amount of green STS probe: 79-99ng/test
Amount of green DXZ1 probe: 8-10ng/test
The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

Counterstain:

150µl per vial (15 tests)
The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius[®] kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all fluorophores and DAPI simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative that should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

FISH Protocol

(Note: Please ensure that exposure of the probe to laboratory lights is limited at all times).

Slide preparation

1. Spot the cell sample onto a glass microscope slide. Allow to dry.
2. Immerse the slide in 2xSSC for 2 minutes at room temperature (RT) without agitation.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 minutes at RT.
4. Allow to dry.

Pre-Denaturation

5. Remove the probe from the freezer and allow it to warm to RT.
6. Ensure that the probe solution is uniformly mixed with a pipette.
7. Remove 10µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
8. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
9. Spot 10µl of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

10. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2 minutes.

Hybridisation

11. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

12. Remove the coverslip and all traces of glue carefully.
13. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
14. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
15. Drain the slide and apply 10µl of DAPI antifade onto each sample.
16. Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
17. View with a fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at/or below RT.

Procedural Recommendations

1. Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
2. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCELL Ltd.
3. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators as these temperatures are critical for optimum product performance.
4. The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
5. Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

Expected Results

In a normal male (XY) cell, there should be one red and two green signals (1R, 2G) or one yellow fusion and one green signal (1Y, 1G). A cell with a deletion of KAL1 and STS regions should have one control green signal (0R, 1G), a cell with STS deletion should have one red and one green signal (1R, 1G) and a cell with KAL1 deletion should have two green signals (0R, 2G).

In a normal female (XX) cell, there should be two red and four green (2R, 4G) or two yellow and two green signals (2Y, 2G). A cell with a deletion of KAL1 and STS regions should have one yellow signal and two green controls (1Y, 2G), a cell with STS deletion should have one yellow, one red and two green signals (1Y, 1R, 2G) and a cell with KAL1 deletion should have one yellow and three green signals (1Y, 3G).

The STS and KAL1 region probes may show strong cross-hybridisation to the Y chromosome, because of sequence homologies in this region.

Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

Additional Information

For additional product information please contact the CytoCELL Technical Support Department.

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.
Cytocell Ltd.