A Sysmex Group Company

## Instructions For Use

REF: LPS 100

### Aquarius® Tissue Pretreatment Kit

PROFESSIONAL USE ONLY

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Further information available at [www.cytocell.com](http://www.cytocell.com)

For use in heat pretreatment and enzyme digestion of Formalin-Fixed, Paraffin-Embedded (FFPE) tissue prior to Fluorescence In Situ Hybridization (FISH) or Chromogenic In Situ Hybridization (CISH) detection.

#### Intended Use

For In Vitro Diagnostic Use. CAUTION: Not for human or animal therapeutic use. Uses other than the labeled intended use may be a violation of local law.

#### Materials Provided

**Reagent 1:** One Litre of Heat Pretreatment Solution, pH 7.0 (Ready to Use)

**Reagent 2:** One 10ml bottle of Enzyme Reagent (Ready to Use)

#### Storage

Store at 2-8°C.

#### Stability

There are no obvious signs to indicate instability of these reagents. They have been quality controlled to assure consistent and reliable performance. Do not use after the expiration date stamped on container. With proper storage there is no significant loss of performance. If the reagents are stored under any conditions other than those specified, conditions must be validated by the user.

#### Tissue Preparation

Slides should be treated with an adhesive before mounting tissue section. Deparaffinize slides and rehydrate tissue sections. Perform heat pretreatment and enzyme digestion as recommended in the instructions that come with the ISH probe. The following procedure may be used if probe instructions are unavailable.

#### Heat Pretreatment

- Heat 500ml Tissue Pretreatment Solution (Reagent 1) in a beaker on a hot plate until it is either boiling or 98 - 100°C. Boil slides for 15 minutes (Note: different incubation times may be required depending on tissue fixation. A 15 minute incubation is a recommended starting point).
- Wash in PBS or dH<sub>2</sub>O at room temperature (RT) for 2 x 3 minutes.

#### Enzyme Digestion

- Cover tissue with 100-200µl of Enzyme Reagent (Reagent 2) for 10 minutes at RT.
- Note: depending on tissue fixative used, different incubation times may be required. Excessive digestion will cause loss of nuclei and chromosome structure. Please refer to Troubleshooting section for details.
- Wash in PBS or dH<sub>2</sub>O at RT for 3 x 2 minutes.
- Dehydrate slides in a series of 70%, 85%, 95% and 100% ethanol for 2 minutes each at room temperature, air dry, and proceed to denaturation and hybridization.

#### Troubleshooting

- Throughout the entire procedure, unless otherwise indicated, it is important that the tissue section does not dehydrate.
- Heat Pretreatment (The most critical step for successful performance):** The specimen must be boiled or heated above 98°C for 15 minutes in Heat Pretreatment Solution.
- Enzyme Digestion (A critical step for successful performance):** Different enzyme incubation times (5 - 15 minutes) may be required, depending on tissue type and fixation method. **For most breast tissues, 10 minutes enzyme digestion at RT will produce the best results. Be sure to pre-warm the Enzyme Pretreatment Reagent to RT prior to adding to the tissue section.** Enzyme pretreatment of the specimen should be evaluated immediately at the completion of the protocol. If nuclei are not counterstained and there is an absent or very weak signal, this may be due to nuclear loss as the result of excessive digestion. If nuclei are strongly counterstained but a

signal is absent in the nuclei, this may be due to under-digestion during the pepsin pretreatment. As an alternative, enzyme pretreatment may also be performed at 37°C for 3 - 10 minutes if optimal results are not attained.

- Probe denaturation at a temperature lower than recommended by the protocol may result in a weak or absent signal.
- Hybridization performed for shorter time periods, or stringent washes performed at higher temperatures, than recommended by the protocol may produce a decrease in or complete loss of the signal.

#### Additional Information

For additional product information please contact the Cytocell Technical Support Department.

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W: www.cytocell.com

#### References/Bibliographie/Bibliografia/Literatur/Bibliografia

- Wilkinson, D.G., *In Situ Hybridization, a Practical Approach*. 2nd ed., Oxford university Press, Oxford, (1998).
- Polak, J.M. and McGee, J. *In Situ Hybridization: Principles and Practice*. Oxford University Press, Oxford, UK, (1998).
- Verma, R.S and Babu, A. *Human Chromosomes: Principles and Techniques*. 2nd ed., Health Professions Division, New York, (1995).
- Leitch, A.R. et al. *In Situ Hybridization-A Practical Guide: Royal Microscopy Society Microscopy Handbooks*. Vol 27, Bios Scientific Publishers, Oxford, UK, (1994).

<b>REF</b>	<b>EN:</b> Catalogue number <b>DE:</b> Bestellnummer <b>FR:</b> Référence du catalogue <b>IT:</b> Riferimento di Catalogo <b>ES:</b> Número de catálogo
<b>IVD</b>	<b>EN:</b> <i>In vitro</i> diagnostic device <b>DE:</b> <i>In-vitro</i> -Diagnostikum <b>FR:</b> Dispositif médical de diagnostic <i>in vitro</i> <b>IT:</b> Dispositivo medico-diagnostico <i>in vitro</i> <b>ES:</b> Producto sanitario para diagnóstico <i>in vitro</i>
<b>LOT</b>	<b>EN:</b> Batch code <b>DE:</b> Loscode <b>FR:</b> Code du lot <b>IT:</b> Codice di lotto <b>ES:</b> Código
	<b>EN:</b> Consult instructions for use <b>DE:</b> Gebrauchsanweisung beachten <b>FR:</b> Consulter la notice d'utilisation <b>IT:</b> Consultare le istruzioni per l'uso <b>ES:</b> Consúltense las instrucciones de uso
	<b>EN:</b> Manufacturer <b>DE:</b> Hersteller <b>FR:</b> Fabricant <b>IT:</b> Fabricante <b>ES:</b> Fabricante
	<b>EN:</b> Use by <b>DE:</b> Verwendbar bis <b>FR:</b> Utiliser jusqu'au <b>IT:</b> Utilizzare entro <b>ES:</b> Fecha de caducidad
	<b>EN:</b> Temperature limitation <b>DE:</b> Temperaturbegrenzung <b>FR:</b> Limites de température <b>IT:</b> Limiti di temperatura <b>ES:</b> Limitación de temperatura
	<b>EN:</b> Sufficient for <n> tests <b>DE:</b> Ausreichend für <b>FR:</b> Suffisant pour <b>IT:</b> Sufficiente per <b>ES:</b> Válido para
<b>CONT</b>	<b>EN:</b> Contents <b>DE:</b> Inhalt <b>FR:</b> Contenu <b>IT:</b> Contenuto <b>ES:</b> Contenido

#### Patents and Trademarks

Aquarius and Cytocell are registered trademarks of Cytocell Ltd.

#### Cytocell Ltd.

Oxford Gene Technology,  
418 Cambridge Science Park,  
Milton Road,  
Cambridge, CB4 0PZ, UK  
T: +44(0)1223 294048  
F: +44(0)1223 294986  
E: probes@cytocell.com  
W: www.cytocell.com

#### LABEL LICENSE

Purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product only for life science research or human diagnostics, as set forth in the accompanying product literature. The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information or data to an unaffiliated third party for payment other than human diagnostics testing services; (3) for therapeutic or prophylactic purposes; (4) to resell, sell or otherwise transfer this product or its components to any third party, or use for any use other than use in life science research, whether or not such product or its components are resold for use in research or human diagnostics.