Optimised, 1-day hybridisation-based NGS protocol yields 1% variant detection in MPN samples, as quickly and cost-effectively as multiplex PCR enrichment

- type]), JAK2 exon 12, MPL W515K/L and CALR exon 9 insertion and deletions.
- [VAF]) of somatic variants of these key MPN driver mutations.
- kit (LPK).

['] Preparation of purified DNA to sequencer-ready libraries in 7 hours 45 minutes

- offers a streamlined, 1-day workflow from purified DNA sample to sequencer.
- poor uniformity of coverage.



Figure 1: Workflow of SureSeq NGS library preparation, from DNA to sequencer.

Panel content

exon 12), MPL (W515K/L/R/A and S505N) and CALR (exon 9).

Study design

- containing variants for each of the targeted regions. Sequencing was conducted on a MiSeq® using a V2 300 bp cartridge (Illumina).

Oxford Gene Technology, Begbroke Science Park, Woodstock Road, Begbroke, Oxfordshire, OX5 1PF, UK T: +44(0)1865 856826 (US: 914-467-5285) E: contact@ogt.com W: www.ogt.com

Graham Speight¹, Ephrem Chin¹, Lyudmila Georgieva¹, Nick Cross² and David Cook¹ ¹Oxford Gene Technology, Oxford, UK and ² National Genetics Reference Laboratory - Wessex, UK

Wild type reads	Variant reads	
1038	157	
1623	78	
2172	23	
2414	9	
2555	1	

Figure 2: Detection of a 52 bp deletion (exon 9 CALR). Wild-type sample (top panel) is compared to a 52 bp somatic deletion (bottom panel).

Figure 3: Detection of a 5 bp deletion (exon 12 JAK2). Wild-type sample (top panel) is compared to a 5 bp somatic deletion (bottom panel).

Technical validation of the Core MPN Panel and enhanced library preparation approach with 14 research samples

#	Gene	Variant	Mean target coverage	% VAF
1	CALR	K385NfsTer47 (5 bp insertion in exon 9 <i>CALR</i>)	3640	37.3
2	JAK2	V617F	1263	10.7
3	JAK2	V617F	1352	9.8
4		Control - no variant detected	2655	
5	MPL	W515R	2750	18.1
6	JAK2	V617F	2971	2.7
7	JAK2	V536F537insFHKIRNEDLIL (33 bp insertion in exon 12 <i>JAK2</i>)	1574	3.1
8	CALR	K385NfsTer47	2321	42.8
9	JAK2	N542E543del	1749	29.9
10	CALR	L367WfsTer46 (52 bp deletion exon 9 <i>CALR</i>)	947	16.3
11	MPL	W515L	1421	60.7
12	MPL	W515L	903	3.8
13	CALR	K385NfsTer47 (5 bp ins in exon 9 <i>CALR</i>)	1265	14.0
14		Control - no variant detected	3234	

Table 2: Data generated using the SureSeq Core MPN Panel in combination with the enhanced LPK was 100% concordant with independent findings (National Genetics Reference Laboratory – Wessex, UK).

- VAF.
- 52 bp deletions and 5 bp insertions).
- cost effective and timely manner.

SureSeq: For Research Use Only; Not for Diagnostic Procedures. This document and its contents are © Oxford Gene Technology IP Limited – 2016. All rights reserved. OGT™ and SureSeq™ are trademarks of Oxford Gene Technology IP Limited. The SureSeq NGS Library Preparation Kit was jointly developed between Oxford Gene Technology and Bioline Reagents Limited.



Results

• Data presented here are from 14 research samples (including 2 wild-type controls) that were processed using the enhanced LPK in combination with the Core MPN Panel.

• See below (Table 2) for the range of variants and VAFs detected from 14 research samples. These include SNVs as well as 5 bp insertions in JAK2 (exon 12) and deletions of up to 52 bp in CALR (exon 9) (Table 2). No variants were identified in the control samples.

Conclusions

• We have successfully utilised the OGT 1-day hybridisation-based SureSeq LPK protocol in combination with the SureSeq Core MPN Panel to reliably and routinely detect somatic SNVs by NGS down to a 1%

• The uniformity of coverage of this approach permitted the detection of CALR and JAK2 indels (including

• This enhanced protocol incorporates an enzymatic fragmentation which permits the high-throughput preparation of 48 samples from genomic DNA to sequencer in a 1-day workflow.

• To achieve >1000x de-duplicated depth (required for confident detection of 1% VAF), 48 samples can be reliably sequenced in a single MiSeq (V2 300 bp) run. This allows the generation of high quality data in a