A Sysmex Group Company

Instructions For Use

REF:
LPT MRK / LPT MRK-S,
LPT MRKxx / LPT MRKxx-S

TeloMark

PROFESSIONAL USE ONLY

ENGLISH/FRANÇAIS/ITALIANO/DEUTSCH/ESPAÑOL
Further information available at www.ogt.com

Intended Use

CytoCell's TeloMark kit is intended to identify, via Fluorescence *In Situ* Hybridisation (FISH), copy number changes or rearrangements involving the subtelomere regions of human chromosomes in patients with mental retardation or developmental delay. Centromere and locus specific probes are intended for enumeration only. It is intended to be used on cultured peripheral blood lymphocytes. This product is for professional use only and is intended to be an adjunct to Clinical Cytogenetics.

Background

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Information

Subtelomere FISH analysis is a commonly used as an adjunct to routine cytogenetic testing in order to detect small rearrangements involving the telomeric regions of chromosomes. Patients can be referred for telomere analyses for a number of different reasons, including genetic disease, autistic disorders, unexplained mental retardation/developmental delay, recurrent miscarriages and haematological malignancies¹⁻⁴.

CytoCell's TeloMark kit consists of 41 subtelomere specific probes, three centromere and six locus specific probes (50 different probes in total). The subtelomere probes represent all chromosome ends apart from the p-arms of the five acrocentric chromosomes and with no distinction between the X and Y due to clones being located in the pseudoautosomal regions. All of these probes map to unique regions of the chromosomes, within 850kb of the true telomere.

The probes in TeloMark are provided as 15 separate mixes, either provided individually or as a kit containing all 15. All probes are directly labelled in orange, green, yellow (orange and green combined) or blue.

Probe Specification

| Mix | Probes | Concentrations |
|--------|---|--|
| Mix 01 | 1p in Green 1q in Orange Xp/Yp in Orange and Green (Yellow) DXZ1 in Blue | 41-62ng/test 18-22ng/test 118-168ng/test 25-37ng/test |
| Mix 02 | 2p in Green 2q in Orange Xq/Yq in Orange and Green (Yellow) DXZ1 in Blue | 41-62ng/test 18-22ng/test 59-84ng/test 25-37ng/test |
| Mix 03 | 3p in Green 3q in Orange 22q in Orange and Green (Yellow) BCR (22q11) in Blue | 41-62ng/test 18-22ng/test 59-84ng/test 76-110ng/test |
| Mix 04 | 4p in Green 4q in Orange 21q in Orange and Green (Yellow) AML1 (RUNX1) (21q22) in Blue | 41-62ng/test 18-22ng/test 59-84ng/test 76-110ng/test |
| Mix 05 | 5p in Green 5q in Orange | 41-62ng/test 18-22ng/test |

| | | |
|--------|--|---|
| Mix 06 | 6p in Green 6q in Orange 13q in Orange and Green (Yellow) 13q14 in Blue | 41-62ng/test 18-22ng/test 59-84ng/test 102-146ng/test |
| Mix 07 | 7p in Green 7q in Orange 14q in Orange and Green (Yellow) TCRAD (14q11.2) in Blue | 41-62ng/test 18-22ng/test 59-84ng/test 76-110ng/test |
| Mix 08 | 8p in Green 8q in Orange 17p in Orange and Green (Yellow) D17Z1 in Blue | 41-62ng/test 18-22ng/test 79-111ng/test 25-37ng/test |
| Mix 09 | 9p in Green 9q in Orange 17q in Orange and Green (Yellow) D17Z1 in Blue | 41-62ng/test 18-22ng/test 59-84ng/test 25-37ng/test |
| Mix 10 | 10p in Green 10q in Orange 15q in Orange and Green (Yellow) PML (15q24) in Blue | 41-62ng/test 18-22ng/test 118-168ng/test 76-110ng/test |
| Mix 11 | 11p in Green 11q in Orange 18p in Orange and Green (Yellow) D18Z1 in Blue | 27-41ng/test 12-14ng/test 45-63ng/test 25-37ng/test |
| Mix 12 | 12p in Green 12q in Orange 18q in Orange and Green (Yellow) D18Z1 in Blue | 82-123ng/test 36-43ng/test 99-139ng/test 25-37ng/test |
| Mix 13 | 16p in Green 16q in Orange | 41-62ng/test 18-22ng/test |
| Mix 14 | 19p in Green 19q in Orange E2A (TCF3) (19p13) in Blue | 41-62ng/test 19-23ng/test 36-51ng/test |
| Mix 15 | 20p in Green 20q in Orange | 21-31ng/test 13-15ng/test |

Materials Provided

Probe: 15µl per vial (5 tests) or 30µl per vial (10 tests)

The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

Counterstain:

150µl per tube (LPT MRKxx/LPT MRKxx-S), 300µl per tube (LPT MRK-S) or 600µl per tube (LPT MRK)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.
6. Take care not to let coverslips touch on the slide after the addition of the probe as this may cause cross contamination between probes.
7. Users of this product must be capable of visually distinguishing between the colours red, blue and green.

Storage and Handling

The kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench top centrifuge.
10. Microscope slides.
11. 12mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.

Fluorescence Microscope Recommendation


For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filters DAPI/FITC/Texas Red or DAPI/FITC/TRITC are optimal for viewing all fluorophores and DAPI simultaneously. Alternatively for viewing orange and green fluorophores use dual bandpass filters FITC/Texas Red or FITC/TRITC. The blue fluorophore has specificity to the Aqua and DEAC spectrum (single bandpass Aqua or DEAC filter is required).

The fluorescence microscope should be checked before use to ensure it is operating correctly. Immersion oil should be suitable for use in fluorescence microscopy and formulated for low autofluorescence. Manufacturers' recommendations should be followed in regards to the life of the lamp and the age of the filters.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative that should be prepared according to the laboratory or institution guidelines.

1. Spot the cell sample onto 3 glass microscope slides with 5 target areas per slide in order to hybridise all 15 TeloMark mixes. Apply 3µl of cell sample per area. Allow to dry.

Example of slide number 1 containing targets 1, 2, 3, 4 and 5. Repeat the same process for slides 2 and 3 with 5 targets per slide.

Recommended slide Pretreatment:

1. Immerse the slide in 2xSSC for 2 minutes at 37°C.
2. Place the slide in freshly made pepsin working solution (5mg of pepsin added to 100ml of 0.01M HCl) for 10 minutes at 37°C.
3. Immerse the slide in phosphate buffered saline (PBS) at RT for 5 minutes.
4. Immerse the slide in post fixation solution (0.95% formaldehyde: 1.0ml of 37% formaldehyde, 0.18g of MgCl₂ and 39.0ml of PBS) for 5 minutes at RT.
5. Immerse the slide in PBS at RT for 5 minutes.
6. Immerse the slide in 70% ethanol at RT. Allow the slide to stand in the ethanol wash for 1 minute.
7. Remove the slide from 70% ethanol. Repeat step 6 with 85% ethanol, followed by 100% ethanol.
8. Allow to air dry.

FISH Protocol

(Note: Please ensure that exposure of the probe to laboratory lights is limited at all times).

Slide preparation (skip this step if the slide was pretreated according to the protocol above)

1. Immerse the slide in 2xSSC for 2 minutes at room temperature (RT) without agitation.
2. Dehydrate in an ethanol series (70%, 85% and 100%) each for 2 minutes at RT.
3. Allow to dry.

Pre-Denaturation

1. Remove the probe from the freezer and allow it to warm to room temperature (RT). Briefly centrifuge tubes before use.
2. Ensure that the probe solution is uniformly mixed with a pipette.
3. Place the probe mixes and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
4. Spot 3µl of probe mixture onto the corresponding target area and carefully apply a 12mm coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

5. Denature the sample and probe simultaneously by heating the slides on a hotplate at 75°C (+/- 1°C) for 2 minutes.

Hybridisation

6. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

7. Remove the coverslip and all traces of glue carefully.
8. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
9. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
10. Drain the slide and apply 10µl of DAPI antifade onto either end of each slide.
11. Cover with a 22x50mm coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
12. View with a fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at/or below RT.

Procedural Recommendations

1. Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
2. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCell Ltd.
3. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators as these temperatures are critical for optimum product performance.
4. The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
5. Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

Expected Results

Normal cells should show the results listed in the table below. For mixes 1 and 2, the number of normal blue signals will vary between male (1 blue) and female (2 blue) samples.

Any deviations from the results listed below should be classed as abnormal.

| Mix | Normal Expected Results |
|--------|---|
| Mix 01 | 2 Green, 2 Orange, 2 Yellow, 1 or 2 Blue (2G, 2O, 2Y, 1/2B) |
| Mix 02 | 2 Green, 2 Orange, 2 Yellow, 1 or 2 Blue (2G, 2O, 2Y, 1/2B) |
| Mix 03 | 2 Green, 2 Orange, 2 Yellow, 2 Blue (2G, 2O, 2Y, 2B) |
| Mix 04 | 2 Green, 2 Orange, 2 Yellow, 2 Blue (2G, 2O, 2Y, 2B) |
| Mix 05 | 2 Green, 2 Orange (2G, 2O) |
| Mix 06 | 2 Green, 2 Orange, 2 Yellow, 2 Blue (2G, 2O, 2Y, 2B) |
| Mix 07 | 2 Green, 2 Orange, 2 Yellow, 2 Blue (2G, 2O, 2Y, 2B) |
| Mix 08 | 2 Green, 2 Orange, 2 Yellow, 2 Blue (2G, 2O, 2Y, 2B) |
| Mix 09 | 2 Green, 2 Orange, 2 Yellow, 2 Blue (2G, 2O, 2Y, 2B) |
| Mix 10 | 2 Green, 2 Orange, 2 Yellow, 2 Blue (2G, 2O, 2Y, 2B) |
| Mix 11 | 2 Green, 2 Orange, 2 Yellow, 2 Blue (2G, 2O, 2Y, 2B) |

| | |
|--------|--|
| Mix 12 | 2 Green, 2 Orange, 2 Yellow, 2 Blue (2G, 2O, 2Y, 2B) |
| Mix 13 | 2 Green, 2 Orange (2G, 2O) |
| Mix 14 | 2 Green, 2 Orange, 2 Blue (2G, 2O, 2B) |
| Mix 15 | 2 Green, 2 Orange (2G, 2O) |

Known Cross Reactivity

| Mix | Known Cross Reactivity |
|--------|-----------------------------------|
| Mix 01 | None |
| Mix 02 | None |
| Mix 03 | 22q with 2q (pericentromeric) |
| Mix 04 | None |
| Mix 05 | None |
| Mix 06 | None |
| Mix 07 | 14q with 16 (centromeric) |
| Mix 08 | 8p to XpYp telomere region |
| Mix 09 | 9q to chromosome 1 (interstitial) |
| Mix 10 | None |
| Mix 11 | 11p to 17p |
| Mix 12 | None |
| Mix 13 | None |
| Mix 14 | None |
| Mix 15 | 20q to 6p and 4q (interstitial) |

Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone. Failure to adhere to the protocol may alter the performance of the assay and may yield erroneous results.

Additional Information

For additional product information please contact the CytoCell Technical Support Department.
T: +44 (0)1223 294048
E: techsupport@cytozell.com
W: www.ogt.com

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Patents and Trademarks
CytoCell is a registered trademark of CytoCell Ltd.
This product contains technology licensed from Life Technologies Corporation