The use of a hybridisation-based NGS enrichment panel for the confident identification of a broad range of low-frequency variants from as little as 50 ng of challenging clinical research FFPE samples

Graham Speight¹, Ephrem Chin¹, Jacqueline Chan¹, Robert Zeillinger², Nicole Concin³, David Cook¹ ¹Oxford Gene Technology, Oxford, UK; ²Medical University of Vienna, Austria and ³Medical University, Dept. of Gynaecology and Obstetrics, Innsbruck, Austria

Introduction

One of the challenges in the treatment of cancer is the high level of genetic complexity and tumour heterogeneity. Detailed information about the genetic profile of each individual tumour can help guide treatment strategies¹. Epithelial ovarian cancer (EOC) is the most lethal gynaecological malignancy², and the type II tumours accounting for approximately 75% of all EOCs, are nearly always detected in advanced stages. These highly aggressive tumours are characterised by their morphological and molecular homogeneity and often (>80% of cases) contain TP53 mutations.

The GANNET53 (Ganetespib in metastatic, p53 mutant, platinum-resistant ovarian cancer) trial started in October 2015 and aims to improve the prognosis and quality of life in platinum-resistant EOC patients (www.gannet53.eu). This European multi-centre clinical trial is currently in stage II during which biomaterials have been collected and analysed using a SureSeqTM hybridisation-based enrichment panel for targeted next-generation sequencing (NGS), to determine the TP53 mutation status of the samples.

Formalin damage in DNA can be reduced using the SureSeq FFPE DNA Repair Mix

Genetic profiling of solid tumours is often problematic as tissue biopsies are typically archived as formalin-fixed, paraffin-embedded (FFPE) blocks, which preserve tissue morphology and permits long-term storage at room temperature. However, the methods used for fixation significantly damage and compromise the quality of nucleic acids in these samples. Formalin damage can fragment the DNA which can impair PCR (a required process in NGS library preparation). Consequently, library yields and the ability to obtain high-quality, meaningful sequence data are compromised, affecting the confident identification of variants.

We tested a range of FFPE derived DNA and found treatment with the SureSeq FFPE DNA Repair Mix significantly improved the mean target coverage, thereby increasing the confidence of the variant calls (Figure 1A). Use of the Repair Mix also enabled a reduction in the amount of DNA input down to 50 ng whilst maintaining good depth of coverage (Figure 1B).

Use of the SureSeq Ovarian Cancer Panel with Type II EOC samples identifies variants in other DNA repair genes in addition to TP53

All 32 samples presented here were provided by the GANNET53 trial and contained ≥40% tumour content. DNA was extracted from tissue curls using standard methods. The resultant DNA was analysed using the hybridisation-based SureSeq Ovarian Cancer Panel and subsequently sequenced on an Illumina MiSeq®. Each sample received between 1/10th to 1/16th of a lane. The average depth of coverage (after removal of PCR duplicates) over the seven target genes (ATM, ATR, BRCA1, minor allele frequencies (MAF) ranging from 1.1 - 79.6%.

In addition to the mutations in TP53, 22 of the samples were found to have additional variants in BRCA1 and BRCA2, (Integrated Genome Viewer [IGV⁴] images of examples are shown in Figures 2 & 3), of which 8 were likely germline (defined as having a MAF between 45-55% or >95%). The remaining putative somatic variants had MAFs ranging from 2.3 to 71.3%. Variants were also found in ATM in 5 of the mutant *TP53* tumour samples (example shown in Figure 4).



Oxford Gene Technology, Begbroke Science Park, Woodstock Road, Begbroke, Oxfordshire, OX5 1PF, UK T: +44(0)1865 856826 (US: 914-467-5285) E: contact@ogt.com W: www.ogt.com

Enrichment methods

Multiplex PCR-based approaches are often used to analyse low input DNA from FFPE samples. However, they are not able to (without molecular barcodes) fully elucidate the true allele frequency and complexity, which is essential to fully evaluate highly heterogenic tumour samples. Hybridisation-based enrichment approaches typically demonstrate better uniformity, are more likely

to preserve the complexity of the original sample, are more tolerant of both DNA quality and unknown variants in the capture region.³

 Table 1: Performance

 comparison of ampliconand hybridisation-based capture methods

Multiplex PCRbased enrichment

Hybridisationbased enrichment



Figure 1: Example data obtained using FFPE DNA extracted from Ovarian tumour samples. Panel A shows that the SureSeq FFPE DNA Repair Mix improves on-target reads (OTR); Panel B demonstrates the use of lower DNA inputs whilst maintaining depth of coverage.

SureSeq: For Research Use Only; Not for Diagnostic Procedures. This document and its contents are © Oxford Gene Technology IP Limited – 2016. All rights reserved. OGT™ and SureSeq™ are trademarks of Oxford Gene Technology IP Limited. The SureSeq NGS Library Preparation Kit was jointly developed between Oxford Gene Technology and Bioline Reagents Limited.

Advantages	Disadvantages
Simple and fast workflow – <1 day Permits low DNA inputs Low start-up costs	 Cannot remove PCR duplicates/bias – obscures true complexity Allelic drop-out due to variants in priming sites Poor uniformity of coverage
Highly uniform coverage Tolerant of variants throughout target region High sensitivity	 Requires greater DNA input Multi-step workflow – 1-2 days



Hybridisation-based enrichment workflow

The SureSeq hybridisation-based workflow was used to determine the genetic profile of 32 ovarian tumours. The workflow of this approach is outlined in Figure 5.



Figure 5: OGT SureSeq workflow. The SureSeq Ovarian Cancer panel targets seven key genes – ATM, ATR, BRCA1, BRCA2, NF1, PTEN and TP53

- with the SureSeq Ovarian Panel.
- tumours.
- pathology, 136(4), pp.527-539.
- Kurman, R.J. and Ie-Ming, S., 2010. The origin and pathogenesis of epithelial ovarian cancer: a proposed unifying theory. The American journal of surgical pathology, 34(3), pp.433-443. Samorodnitsky, E., et al., 2015. Evaluation of Hybridization Capture Versus Amplicon-Based Methods for Whole-Exome
- Sequencing. *Human mutation*, 36(9), pp.903-914.



Conclusions

• We have demonstrated that it is possible to gain useful genetic information from as little as 50 ng of DNA derived from FFPE material when using the SureSeq FFPE DNA Repair Mix in combination

• Using this hybridisation-based capture approach we were able to identify a range of germline and somatic mutations in TP53 from ovarian tumour tissue collected as part of the GANNET53 trial. Some of these mutations are likely to be the oncogenic drivers in the pathogenesis of these

• We also identified numerous BRCA1 and BRCA2 mutations in 71% of the samples confirming independent observations that these genes along with ATM are important in ovarian carcinomas.

References

Ross, J.S. and Cronin, M., 2011. Whole cancer genome sequencing by next-generation methods. American journal of clinical

Thorvaldsdóttir, H., Robinson, J.T. and Mesirov, J.P., 2013. Integrative Genomics Viewer (IGV): high-performance genomics data visualization and exploration. *Briefings in bioinformatics*, 14(2), pp.178-192.