A Sysmex Group Company

Instructions For Use
REF: LPH 087-S / LPH 087

CLL Plus Screening Panel

PROFESSIONAL USE ONLY

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Further information available at www.cytocell.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Information

Alpha Satellite 12 Plus in Red

Trisomy of chromosome 12 is a relatively common chromosome abnormality in B-CLL, being present in around 16% of patients, when evaluated by FISH on interphase and metaphase cells¹ and is now thought to have strong prognostic significance².

13q14.3 Deletion Probe

Deletions of 13q14 are a common event in CLL. This region is found to be heterozygously deleted in 30-60% and homozygously deleted in 10-20% of CLL patients³.

P53 (TP53) Deletion Probe

P53 (TP53) is a tumour suppressor gene that mediates the apoptosis of damaged cells. It has been found that approximately 17% of B-CLL patients have deletions of the P53 gene⁴. Patients exhibiting this genotype are associated with a poor prognosis⁵ as they harbour a proliferating population of damaged cells.

ATM Deletion Probe

The ATM gene is an important regulatory candidate in cell damage management. When it is deleted, damaged cells are neither repaired nor apoptosed and are allowed to proliferate. Screening for deletions of ATM is important to allow informed therapy choices for CLL patients, especially as deletions confer poor prognosis⁶.


MYB Deletion probe

Deletions of chromosome 6q are found in B-CLL, with breakpoints being reported in bands 6q13, q15 or q21. The MYB gene is essential in haematopoietic cell proliferation and differentiation. It is located in band 6q23.3 and is provided as a marker for 6q deletion as it is distal to reported deletion breakpoints (6q13, 6q15 and 6q21)⁶.

Probe Specification

13q14.3 Deletion Probe


13q14.2-q14.3, Red
13qter, 13q34, Green

The 13q14.2-q14.3 probe, labelled in red, covers the D13S319 and D13S25 markers. The 13qter subtelomere specific probe (clone 163C9), labelled in green, allows identification of chromosome 13 and acts as a control probe.

ATM Deletion Probe

ATM, 11q22.3, Red
D11Z1, 11p11.1-q11.1, Green

The ATM probe is 182kb, labelled in red, and covers the telomeric end of the NPAT gene and the centromeric end of the ATM gene to just beyond the D11S3347 marker. The probe mix also contains a control probe for the 11 centromere (D11Z1) labelled in green.

Alpha Satellite 12 Plus Probe in Red


D12Z3, 12p11.1-q11.1, Red

The Alpha Satellite 12 Plus Probe is a repeat sequence probe, labelled in red, which recognizes the centromeric repeat sequence D12Z3.


P53 Deletion Probe

P53, 17p13.1, Red
D17Z1, 17p11.1-q11.1, Green

The P53 probe is 161kb, labelled in red, covers the whole P53 (TP53) gene, extending 66kb telomeric to the gene and covering a region centromeric to the gene, to just beyond the marker D17S655. The probe mix also contains a control probe for the 17 centromere (D17Z1) labelled in green.

MYB Deletion Probe
MYB, 6q23.3, Red
D6Z1, 6p11.1-q11.1, Green

The MYB probe is 183kb, labelled in red, covers the entire MYB gene and a region telomeric to the gene, 137kb beyond the marker AFMA074ZG9. The probe mix also contains a control probe for the 6 centromere (D6Z1) labelled in green.

Materials Provided

13q14.3 Deletion Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of red 13q14.3 probe: 24-30ng/test

Amount of green 13qter probe: 72-90ng/test

ATM Deletion Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of red ATM probe: 40-50ng/test

Amount of green D11Z1 probe: 30-37.5ng/test

Alpha Satellite 12 Plus in Red: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of red D12Z3 probe: 5-6ng/test

P53 Deletion Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of red P53 probe: 13-16ng/test

Amount of green D17Z1 probe: 7-9ng/test

MYB Deletion Probe: 50µl per vial (5 tests) or 100µl per vial (10 tests)

Amount of red MYB probe: 15-19ng/test

Amount of green D6Z1 probe: 12-15ng/test

The probes are provided premixed in hybridisation solution (Formamide; Dextran Sulphate; SSC) and are ready to use.

Counterstain: 500µl per vial (50 tests) or 600µl per vial (60 tests)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
4. DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Aquarius® kit should be stored between -25°C to -15 °C in a freezer until the expiry date indicated on the kit label. The probe and counterstain vials must be stored in the dark.

Equipment Necessary but not Supplied

1. Hotplate (with a solid plate and accurate temperature control up to 80°C).
2. Variable volume micropipettes and tips range 1µl - 200µl.
3. Water bath with accurate temperature control at 72°C.
4. Microcentrifuge tubes (0.5ml).
5. Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
6. Plastic or glass coplin jars.
7. Forceps.
8. Fluorescence grade microscope lens immersion oil.
9. Bench top centrifuge.
10. Microscope slides.
11. 24x24mm coverslips.
12. Timer.
13. 37°C incubator.
14. Rubber solution glue.

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all fluorophores and DAPI simultaneously.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells or cultured bone marrow cells fixed in Carnoy's fixative that should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

FISH Protocol

(Note: Please ensure that exposure of the probe to laboratory lights is limited at all times).

Slide preparation

1. Spot the cell sample onto a glass microscope slide. Allow to dry.
2. Immerse the slide in 2xSSC for 2 minutes at room temperature (RT) without agitation.
3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 minutes at RT.
4. Allow to dry.

Pre-Denaturation

5. Remove the probe from the freezer and allow it to warm to RT.
6. Ensure that the probe solution is uniformly mixed with a pipette.
7. Remove 10µl of probe per test, and transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
8. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
9. Spot 10µl of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

10. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2 minutes.

Hybridisation

11. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

12. Remove the coverslip and all traces of glue carefully.
13. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
14. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
15. Drain the slide and apply 10µl of DAPI antifade onto each sample.
16. Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
17. View with a fluorescence microscope.

Stability of Finished Slides

FISHed slides remain analysable for up to 1 month if stored in the dark at/or below RT.

Procedural Recommendations

1. Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
2. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by Cytocell Ltd.
3. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators as these temperatures are critical for optimum product performance.
4. The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
5. Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

Expected Results

13q14.3 Deletion Probe

In a normal cell there should be two red and two green signals (2R, 2G). A cell with a hemizygous deletion of the 13q14.3 should have one red and two green signals (1R, 2G) whilst a cell with a homozygous deletion should have no red and two green signals (0R, 2G).

ATM Deletion Probe

In a normal cell there should be two red and two green signals (2R, 2G) whilst a cell with an ATM deletion should have one red and two green signals (1R, 2G).

Alpha Satellite 12 Plus in Red

In a normal cell two red signals should be observed (2R). In a cell with trisomy 12 there should be three red signals (3R).

P53 Deletion Probe

In a normal cell there should be two red and two green signals (2R, 2G) whilst a cell with a P53 deletion should have one red and two green signals (1R, 2G).

MYB Deletion Probe

In a normal cell there should be two red and two green signals (2R, 2G) whilst a cell with a MYB deletion should have one red and two green signals (1R, 2G).

Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone.

Additional Information

For additional product information please contact the Cytocell Technical Support Department.

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Patents and Trademarks
Aquarius and Cytocell are registered trademarks of Cytocell Ltd.