

# CytoSure



A Sysmex Group Company

## Comprehensive FH NGS Panel



## Features

### Detection of CNVs as well as SNVs with a single assay

- Enabled by the exon resolution of the targeted genes

### Pre-optimised content that meets your technical requirements

- No more laborious in-house optimisation, decreasing assay development time

### Bespoke panel content

- Sequence only what's relevant for your research

### Interpret, OGT's complimentary analysis software

- Designed to give unparalleled CNV and SNV calling

### Introduction

Familial Hypercholesterolaemia (FH) is a genetic condition which results in a high cholesterol level and subsequently leads to a higher risk of early heart disease. It affects approximately 1 in 250 people with around 34 million cases worldwide<sup>1</sup>.

OGT is offering an optimised NGS panel which has selected the most relevant genes and SNPs implicated in FH, for your research needs. Together with the complimentary Interpret analytical software, the CytoSure® Comprehensive FH NGS platform provides the optimal solution for FH research.

### Genes at exon-level resolution



### SNPs associated with FH

rs10401969 SUGP1 Intron Variant	rs2479409 PCSK9 : 2KB Upstream Variant
rs10455872 LPA Intron Variant	rs2650000
rs11206510 Ch1	rs3757354 MYLIP : 2KB Upstream Variant
rs11220462 ST3GAL4 intron variant	rs3798220 LPA : Missense Variant
rs121909548 SERPINC1 : Missense Variant	rs3846663 HMGCR : Intron Variant
rs12740374 CELSR2 UTR variant	rs4149015 SLCO1B1 : 2KB Upstream Variant
rs1346268 GATM : Intron Variant	rs4149056 SLCO1B1 : Missense Variant
rs1367117 APOB : Missense Variant	rs428785 ADAMTS1 : Missense Variant
rs1501908	rs429358 APOE : Missense Variant
rs1564348 SLC22A1 : Intron Variant	rs4299376 ABCG8 : Intron Variant
rs1719247	rs4693075 COQ2 : Intron Variant
rs17244841 HMGCR : Intron Variant	rs515135
rs1799768 SERPINE1 : 2KB Upstream Variant	rs6025 F5 : Missense Variant
rs1799963 F2 : 3 Prime UTR Variant	rs6102059 LOC102724968 : Intron Variant
rs1800562 HFE : Missense Variant, LOC108783645 : 2KB Upstream Variant	rs629301 CELSR2 : 3 Prime UTR Variant
rs1801131 MTHFR : Missense Variant	rs6511720 LDLR : Intron Variant, LDLR-AS1 : 2KB Upstream Variant
rs1801133 MTHFR : Missense Variant	rs6544713 ABCG8 : Intron Variant
rs2032582 ABCB1 : Missense Variant	rs7412 APOE : Missense Variant
rs2231142 ABCG2 : Missense Variant	rs8017377 NYNRIN : Missense Variant
rs2306283 SLCO1B1 : Missense Variant	

Some genes/SNPs may not be available in your region - contact OGT for more details.

### CNV and SNV detection from a single assay

The hybridisation enrichment methodology, combined with our bait design expertise, allows generation of panels with outstanding completeness and coverage uniformity. Together, this allows the areas of CNV to be easily identified within each sample using our proprietary algorithm — delivering an increased understanding of the sample without an increase in cost or time (Figure 1).




Figure 1: CNV in *LDLR* gene shown using IGV from the Broad Institute **A**: Red bars indicate areas of CNV (data from aCGH), purple bars represent deleted exons (data from NGS): 5 samples are shown, each with at least one area of CNV. There is complete concordance between the aCGH and NGS data. Note the evenness of the NGS coverage (even peak height) across each exon, allowing the areas of CNV to be easily identified. The data from the custom CytoSure aCGH array, confirms the deletions in *LDLR*: **B**: A 2 exon deletion and **C**: deletion of 2 exons and 4 exons, corresponding to samples 3 and 5 in **A** respectively.

The Interpret software has been designed to easily visualise CNVs and SNVs, with an intuitive interface to switch between different sets of results. Interpret also has simple to use protocols and filtering options, to easily target the results of interest.

The CytoSure Comprehensive FH NGS panel has the ability to detect CNVs in whole genes, at exon resolution (Figures 2–3) and can target select SNPs that have been implicated in FH (Figure 4). The CytoSure Comprehensive FH panel can also detect SNVs and Indels within genes, as demonstrated by Figure 5.



Figure 2: Double deletion on the *LDLR* gene, as visualised by Interpret software.




Figure 3: Duplication on the *LDLR* gene, as visualised by Interpret software.




Figure 4: Deletion on the *SLOCO1B1* gene, as visualised by Interpret software.




Figure 5: Missense variant on the *APOB* gene, as visualised by Interpret software.

### Product Specifications

Targeted Genes	8
Targeted SNPs	39
CNV Resolution	Exon Level
SNVs and Indels	Targeted Genes (8)
Mean Target Coverage	> 300x
Recommended DNA Input	> 500 ng high quality DNA
Panel Size	50Kb
Samples per MiSeq® v2 run	16

### Ordering information

UK +44 (0) 1865 856800

US +1 914 467 5285

contact@ogt.com

ogt.com

Product	Contents	Cat. No.
CytoSure Comprehensive FH Panel (16 reactions)	Enrichment baits; Interpret Software	601004-16
SureSeq NGS Library Preparation Complete Solution (16)	Bundle of 1x SureSeq NGS Library Preparation Kit (16), containing adaptors, PCR primers and enzymes, 1x SureSeq NGS Index Kit – Collection A, 1x SureSeq NGS Hyb & Wash Kit (16), 1x Dynabeads M270 Streptavidin (2ml) and 1x AMPure XP beads (10ml). Sufficient for 16 samples	500084
SureSeq NGS Library Preparation Complete Solution (48)	Bundle of 1x SureSeq NGS Library Preparation Kit (48), containing adaptors, PCR primers and enzymes, 1x SureSeq NGS Index Kit – Collection B, 3x SureSeq NGS Hyb & Wash Kit (16), 3x Dynabeads M270 Streptavidin (2ml) and 3x AMPure XP beads (10ml). Sufficient for 48 samples	500085

Talk to us about your requirements, email [contact@ogt.com](mailto:contact@ogt.com) or visit [ogt.com/FH](http://ogt.com/FH).

### References

1. Goldberg AC, and Gidding, SS. Knowing the Prevalence of Familial Hypercholesterolemia Matters. *Circulation*, 2016; 133 (11).



A Sysmex Group Company

**What binds us,  
makes us.**

Oxford Gene Technology Ltd., Begbroke Science Park, Woodstock Road, Begbroke, Oxfordshire, OX5 1PF, UK

**SureSeq™ and CytoSure®:** For Research Use Only; Not for Diagnostic Procedures. This document and its contents are © Oxford Gene Technology IP Limited – 2021. All rights reserved. OGT™, CytoSure® and SureSeq™ are trademarks of Oxford Gene Technology IP Limited. The SureSeq NGS Library Preparation Kit was jointly developed between Oxford Gene Technology and Bionline Reagents Limited. Dynabeads is a trademark of Thermo Fisher Scientific and AMPure® is a registered trademark of Beckman Coulter Inc.