

Guide to expected abnormal clinical results

AML1/ETO (RUNX1/RUNX1T1) Translocation, Dual Fusion FISH Probe Kit

Cat. No. USA LPH-026

Expected normal signal pattern

In a normal cell, two red and two green signals (2R, 2G) are expected.

Expected abnormal signal pattern

In a cell with a t(8;21)(q22;q22) translocation the expected signal pattern will be one red, one green and two fusions (1R, 1G, 2F).

Other signal patterns are possible in aneuploid/unbalanced specimens.

Characterization of Normal Cut-off Values*

Abnormal signal pattern	Number of samples analyzed to generate the cut-off	Number of nuclei evaluated per sample	Maximum number of false positive signal pattern	Normal cut-off value (per 200 nuclei)	Normal cut-off value (%)
1R, 1G, 2F	1290	200	1	5	2.3

*The device has not been specifically validated in patients with <20% blast count.

Del(20q) Deletion FISH Probe Kit

Cat. No. USA LPH-020

Expected normal signal pattern

In a normal cell, two red and two green signals (2R, 2G) are expected.

Expected abnormal signal pattern

One red and one green signal pattern (1R, 1G) will be observed in cells with either monosomy or hemizygous deletion of both bands on 20q.

Other signal patterns are possible in aneuploid/unbalanced specimens.

Characterization of Normal Cut-off Values

Abnormal signal pattern	Number of samples analyzed to generate the cut-off	Number of nuclei evaluated per sample	Maximum number of false positive signal pattern	Normal cut-off value (per 200 nuclei)	Normal cut-off value (%)
1R, 1G	1300	200	6	12	5.7

CBFB (CBFB)/MYH11 Translocation, Dual Fusion FISH Probe Kit

Cat. No. USA LPH-022

Expected normal signal pattern

In a normal cell, two red and two green signals (2R, 2G) are expected.

Expected abnormal signal pattern

In a cell with an inv(16)(p13q22) or a t(16;16)(p13;q22) the expected signal pattern will be one red, one green and two fusions (1R, 1G, 2F).

Other signal patterns are possible in aneuploid/unbalanced specimens.

Characterization of Normal Cut-off Values*

Abnormal signal pattern	Number of samples analyzed to generate the cut-off	Number of nuclei evaluated per sample	Maximum number of false positive signal pattern	Normal cut-off value (per 200 nuclei)	Normal cut-off value (%)
1R, 1G, 2F	1300	200	1	5	2.3

*The device has not been specifically validated in patients with <20% blast count.

EV11 (MECOM) Breakpart FISH Probe Kit

Cat. No. USA LPH-036

Expected normal signal pattern

In a normal cell, two red/green/blue fusion signals (2RGB) are expected.

Expected abnormal signal pattern

In a cell with a t(3;3)(q21;q26.2) translocation, one red signal, one green/blue fusion signal and one red/green/blue fusion signal will be observed (1R, 1GB, 1RGB).

The three-color strategy shows the presence of either a translocation or an inversion and allows for each different type of rearrangement to be distinguished.

Expected normal signal pattern

In a cell with an inv(3)(q21q26.2) inversion, one red/green fusion, one separate blue signal and one red/green/blue fusion signal will be observed (1RG, 1B, 1RGB).

Other signal patterns are possible in aneuploid/unbalanced specimens.

Characterization of Normal Cut-off Values

Abnormal signal pattern	Number of samples analyzed to generate the cut-off	Number of nuclei evaluated per sample	Maximum number of false positive signal pattern	Normal cut-off value (per 200 nuclei)	Normal cut-off value (%)
1R, 1GB, 1RGB	25	200	3	8	4
1RG, 1B, 1RGB	25	200	3	8	4

Del(5q) Deletion FISH Probe Kit

Cat. No. USA LPH-024

Expected normal signal pattern

In a normal cell, two red and two green signals (2R, 2G) are expected.

Expected abnormal signal pattern

A cell with a hemizygous deletion of 5q31.2 will have one red and two green signals (1R, 2G).

Other signal patterns are possible in aneuploid/unbalanced specimens.

Characterization of Normal Cut-off Values

Abnormal signal pattern	Number of samples analyzed to generate the cut-off	Number of nuclei evaluated per sample	Maximum number of false positive signal pattern	Normal cut-off value (per 200 nuclei)	Normal cut-off value (%)
1R, 2G	1300	200	7	13	6.3

MLL (KMT2A) Breakpart FISH Probe Kit

Cat. No. USA LPH-013

Expected normal signal pattern

In a normal cell, two red/green fusion signals are expected (2F).

Expected abnormal signal pattern

In a cell with a balanced MLL (KMT2A) rearrangement, the expected signal pattern will be one red, one green and one fusion (1R, 1G, 1F).

Other signal patterns are possible in aneuploid/unbalanced specimens.

Characterization of Normal Cut-off Values

Abnormal signal pattern	Number of samples analyzed to generate the cut-off	Number of nuclei evaluated per sample	Maximum number of false positive signal pattern	Normal cut-off value (per 200 nuclei)	Normal cut-off value (%)
1R, 1G, 1F	1600	200	3	8	3.8

Del(7q) Deletion FISH Probe Kit

Cat. No. USA LPH-025

Expected normal signal pattern

In a normal cell, two red and two green signals (2R, 2G) are expected.

Expected abnormal signal pattern

One red and one green signal pattern (1R, 1G) will be observed in cells with either monosomy 7 or hemizygous deletion of both CDRs on 7q.

Other signal patterns are possible in aneuploid/unbalanced specimens.

Characterization of Normal Cut-off Values

Abnormal signal pattern	Number of samples analyzed to generate the cut-off	Number of nuclei evaluated per sample	Maximum number of false positive signal pattern	Normal cut-off value (per 200 nuclei)	Normal cut-off value (%)
1R, 1G	1300	200	9	15	7.4

P53 (TP53) Deletion FISH Probe Kit

Cat. No. USA LPH-017

Expected normal signal pattern

In a normal cell, two red and two green signals (2R, 2G) are expected.

Expected abnormal signal pattern

A cell with a TP53 deletion, will have one red and two green signals (1R, 2G).

Other signal patterns are possible in aneuploid/unbalanced specimens.

Characterization of Normal Cut-off Values

Abnormal signal pattern	Number of samples analyzed to generate the cut-off	Number of nuclei evaluated per sample	Maximum number of false positive signal pattern	Normal cut-off value (per 200 nuclei)	Normal cut-off value (%)
1R, 2G	1600	200	8	14	6.8

Further information US +1 833 208 4648 | OGTorders@sysmex.com | ogt.com | Oxford Gene Technology Inc., 520 White Plains Road, Suite 500, Tarrytown, NY 10591, USA

CytoCell®: This document and its contents are © Oxford Gene Technology IP Limited – 2021. All rights reserved. The CytoCell AML/MDS range of FISH probe test kits are fluorescence *in situ* hybridization (FISH) tests used to detect common chromosomal rearrangements in fixed bone marrow specimens from patients with acute myeloid leukemia (AML) or myelodysplastic syndrome (MDS). The tests are indicated for the characterization of patient specimens consistent with World Health Organization guidelines for Classification of Tumours of Haematopoietic and Lymphoid Tissues (Revised 4th Edition) and in conjunction with other clinicopathological criteria. The assay results are to be interpreted by a qualified pathologist or cytogeneticist. The tests are not intended for use as a stand-alone diagnostic, disease screening, or as a companion diagnostic. Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone. Failure to adhere to the protocol may affect the performance and lead to false results. Each lab is responsible for establishing their own cut-off values. Each laboratory should test sufficiently large number of samples to establish normal population distribution of the signal levels and to assign a cut-off value. The product is for professional use only and is intended to be interpreted by a qualified Pathologist or Cytogeneticist. Refer to individual test kit Package Insert for the specific intended use and limitations. For *In Vitro* Diagnostic Use. Rx only. Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representatives for availability. For sale in the US only. This product has not been licensed in accordance with Canadian law.