
Instructions For Use

REF: LPT xxxR/G

Subtelomere Specific Probes

PROFESSIONAL USE ONLY

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Further information available at www.ogt.com

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei from fixed cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Recent developments have meant that this valuable technique can now be applied as an essential diagnostic tool in prenatal, haematological and pathological chromosomal analysis. Target DNA, after fixation and denaturation, is available for annealing to a similarly denatured, fluorescently labelled DNA probe, which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed and the DNA is counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

Probe Information

Chromosomal rearrangements involving the ends of chromosomes have emerged as an important cause of genetic disease given the gene-rich nature of the regions adjacent to the telomeres¹. The importance of such subtelomeric chromosome rearrangements has been clearly shown by their observed association with unexplained mental retardation and congenital abnormalities². Individual subtelomere specific probes have been used to focus on particular subtelomeric regions, and have resulted in the establishment of syndromes such as the chromosome 1p36 deletion syndrome^{3,10} and the 22q13.3 deletion syndrome⁴. The probes are also finding applications in the investigation of autistic disorders⁵, recurrent miscarriages⁶ and haematological malignancies⁷. CytoCell's subtelomere specific probes are located in the most distal region of chromosome specific DNA on each chromosome. Beyond this unique sequence material is the 100 to 300kb region of telomere associated repeat followed by the cap of between 3 to 20kb of tandemly repeated (TTAGGG)_n sequence⁸. The probes have been chosen from the most distal unique sequence to provide the best possible specificity, whilst also being applicable for routine use for the examination of subtelomeric enumeration and integrity. The original second-generation set of probes is derived from PAC clones⁹ and was established in conjunction with the Institute of Molecular Medicine, part of Oxford University, in the UK¹¹. Continuing product improvements have led to some substitutions with alternative cosmid (35-40kb) or BAC (150kb) clones to give improved signal strength or chromosome specificity.

Probe Specification

The subtelomere specific probe range identifies 41 of the 46 human telomeres as it excludes the p-arm telomeres of the acrocentric chromosomes. The X and Y chromosome p-arms share the same subtelomere clone (839D20) as do the X and Y q-arms (C8.2/1 and 225F6) due to the pseudoautosomal nature of these regions. The probes are directly labelled with either a red or a green fluorophore. For detailed probe specifications refer to Table 1.

Table 1: Probe Specifications

Probe	Catalogue Number	Clone name	Marker	Accession Number (if available)
1p	LPT 01PR/G	CEB108	RH120573	-
1q	LPT 01QR/G	160H23	GDB:315525	D1S3739
2p	LPT 02PR/G	dJ892G20	D2S2983	D2S2983
2q NP	LPT 02QNPR/G	172113	D2S447	D2S2986
3p	LPT 03PR/G	dJ1186B18	D3S4559	D3S4559
3q	LPT 03QR/G	196F4	D3S1272	D3S1272
4p	LPT 04PR/G	36P21	D4S3360	D4S3360
4q	LPT 04QR/G	963K6	D4S139	-
5p	LPT 05PR/G	189N21	RH120167	-
5q	LPT 05QR/G	240G13	D5S2907	D5S2907
6p	LPT 06PR/G	62111	STS-H99640	-
6q	LPT 06QR/G	57H24	D6S2522	D6S2522
7p	LPT 07PR/G	109a6	RH104000	RH104000
7q	LPT 07QR/G	2000a5	RH48601	RH48601
8p	LPT 08PR/G	dJ580L5	RH40619	D8S2333
8q	LPT 08QR/G	489D14	D8S595	D8S1925

9p	LPT 09PR/G	43N6	RH65569	RH65569
9q	LPT 09QR/G	112N13	D9S2168	D9S2168
10p	LPT 10PR/G	306F7	STS-N35887	D10S2488
10q	LPT 10QR/G	137E24	RH44494	RH44494
11p	LPT 11PR/G	dJ908H22	D11S2071	D11S2071
11q	LPT 11QR/G	dJ770G7	D11S4974	D11S4974
12p	LPT 12PR/G	496A11	D12S200	D12S200
12q	LPT 12QR/G	221K18	RH81094	D12S2343
13p	LPT 13PR/G	163C9	D13S1825	D13S1825
14q	LPT 14QR/G	dJ820M16	D14S1420	D14S1420
15q	LPT 15QR/G	154P1	D15S936	D15S936
16p	LPT 16PR/G	121I4	SHGC-16929(UCSC)	D16S3400
16q	LPT 16QR/G	240G10	RH80305	RH80305
17p	LPT 17PR/G	2111b1	D17S2199	D17S2199
17q	LPT 17QR/G	362K4	362K4 For and Rev	D17S2200
18p	LPT 18PR/G	74G18	D18S552	D18S552
18q	LPT 18QR/G	dJ964M9	D18S1390	D18S1390
19p	LPT 19PR/G	dJ546C11	D19S676E	-
19q	LPT 19QR/G	F21283	RH102404	RH102404
20p	LPT 20PR/G	dJ1061L1	D20S210	D20S502
20q	LPT 20QR/G	81F12	RH10656	-
21q	LPT 21QR/G	63H24	D21S1446	D21S1575
22q	LPT 22QR/G	99K24	D22S1726	D22S1726
XpYp	LPT XYPR/G	839D20	DXYS129	DXYS129
XqYq	LPT XYQR/G	225F6 C8.2/1	DXYS154 SYBL1	Z43206

R specifies a red label and G specifies a green label

**This probe is specific for the p-arms of both X and Y.

***This probe is specific for q-arms of both X and Y.

This kit contains only one of the probes from the range of directlylabelled subtelomere specific probes.

Materials Provided

Probe: 15µl per vial (5 tests)

Amount of red or green subtelomere specific probe: minimum of 40ng/test

The probe is produced in a concentrated form. It is labelled with either a red or a green fluorophore. The probe is provided in hybridisation solution (Formamide; Dextran Sulphate; SSC).

Hybridisation solution (Formamide; Dextran Sulphate; SSC): 150µl per vial

Counterstain: 150µl per vial (15 tests)

The counterstain is DAPI antifade (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)).

Warnings and Precautions

- For *in vitro* diagnostic use. For professional use only.
- Wear gloves when handling DNA probes and DAPI counterstain.
- Probe mixtures contain formamide, which is a teratogen; do not breathe fumes or allow skin contact. Wear gloves, a lab coat, and handle in a fume hood. Upon disposal, flush with a large volume of water.
- DAPI is a potential carcinogen. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
- All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.
- Users of this product must be capable of visually distinguishing between the colours red, blue and green.

Storage and Handling

Storage and Handling
-25°C -15°C
The probe should be stored between -25°C to -15°C in a freezer until the expiry date indicated on the kit label. The probe vials must be stored in the dark.

Equipment Necessary but not Supplied

- Hotplate (with a solid plate and accurate temperature control up to 80°C).
- Variable volume micropipettes and tips range 1µl - 200µl.
- Water bath with accurate temperature control at 72°C.
- Microcentrifuge tubes (0.5ml).
- Fluorescence microscope (Please see Fluorescence Microscope Recommendation section).
- Plastic or glass coplin jars.
- Forceps.
- Fluorescence grade microscope lens immersion oil.
- Bench top centrifuge.
- Microscope slides.
- 24x24mm coverslips.
- Timer.
- 37°C incubator.
- Rubber solution glue.

Fluorescence Microscope Recommendation

For optimal visualisation of the probe we recommend a 100-watt mercury lamp and plan apochromat objectives x63 or x100. The Triple bandpass filter DAPI/FITC/Texas Red is optimal for viewing all fluorophores and DAPI simultaneously.

The fluorescence microscope should be checked before use to ensure it is operating correctly. Immersion oil should be suitable for use in fluorescence microscopy and formulated for low autofluorescence. Manufacturers' recommendations should be followed in regards to the life of the lamp and the age of the filters.

Sample Preparation

The kit is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative that should be prepared according to the laboratory or institution guidelines. Prepare air dried samples on microscope slides according to standard cytogenetic procedures.

FISH Protocol

(Note: Please ensure that exposure of the probe to laboratory lights is limited at all times).

Slide preparation

- 1. Spot the cell sample onto a glass microscope slide. Allow to dry.
- 2. Immerse the slide in 2xSSC for 2 minutes at room temperature (RT) without agitation.
- 3. Dehydrate in an ethanol series (70%, 85% and 100%), each for 2 minutes at RT.
- 4. Allow to dry.

Pre-Denaturation

- 5. Remove the probe from the freezer and allow it to warm to RT. Briefly centrifuge tubes before use.
- 6. Ensure that the probe solution is uniformly mixed with a pipette.
- 7. Using fresh pipette tips remove (final volume of 10µl of probe solution):
 - for a **single probe hybridisation**: 3µl of probe and 7µl of hybridisation solution per test
 - for a two probe hybridisation: 3µl of each probe and 4µl of hybridisation solution per test
 - for a three probe hybridisation: 3µl of each probe and 1µl of hybridisation solution per testand transfer it to a microcentrifuge tube. Quickly return the remaining probe to the freezer.
- 8. Place the probe and the sample slide to prewarm on a 37°C (+/- 1°C) hotplate for 5 minutes.
- 9. Spot 10µl of probe mixture onto the cell sample and carefully apply a coverslip. Seal with rubber solution glue and allow the glue to dry completely.

Denaturation

- 10. Denature the sample and probe simultaneously by heating the slide on a hotplate at 75°C (+/- 1°C) for 2 minutes.

Hybridisation

- 11. Place the slide in a humid, lightproof container at 37°C (+/- 1°C) overnight.

Post-Hybridisation Washes

- 12. Remove the coverslip and all traces of glue carefully.
- 13. Immerse the slide in 0.4xSSC (pH 7.0) at 72°C (+/- 1°C) for 2 minutes without agitation.
- 14. Drain the slide and immerse it in 2xSSC, 0.05% Tween-20 at RT (pH 7.0) for 30 seconds without agitation.
- 15. Drain the slide and apply 10µl of DAPI antifade onto each sample.
- 16. Cover with a coverslip, remove any bubbles and allow the colour to develop in the dark for 10 minutes.
- 17. View with a fluorescence microscope.

Stability of Finished Slides

FISH slides remain analysable for up to 1 month if stored in the dark at/or below RT.

Procedural Recommendations

- 1. Baking or ageing of slides is not recommended as it may reduce signal fluorescence.
- 2. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by CytoCell Ltd.
- 3. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators as these temperatures are critical for optimum product performance.
- 4. The wash concentrations, pH and temperatures are important as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
- 5. Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

Expected Results for a single probe hybridisation

A normal sample will show a signal in the telomeric region of each of the relevant chromosome homologues, or 2 signals in interphase cells.

Probe	Catalogue Number	Known cross-hybridisations
8p	LPT08PR/G	8p with 1p and 3q
9q	LPT09QR/G	9q with 10p, 16p, 18p and XqYq
11p	LPT11PR/G	11p with 17p
11q	LPT11QR/G	11q with 12q (interstitial)
12p	LPT12PR/G	12p with 6p and 20q
14q	LPT14QR/G	14q with 16 centromere
17q	LPT17QR/G	17q with 1p, 5q, 6q and 11p
19p	LPT19PR/G	19p with 20q
20q	LPT20QR/G	20q with 6p
22q	LPT22QR/G	22q with 2q (interstitial)

Limitations

Reporting and interpretation of FISH results should be consistent with professional standards of practice and should take into consideration other clinical and diagnostic information. This kit is intended as an adjunct to other diagnostic laboratory tests and therapeutic action should not be initiated on the basis of the FISH result alone. Failure to adhere to the protocol may alter the performance of the assay and may yield erroneous results.

Additional Information

For additional product information please contact the CytoCell Technical Support Department.
T: +44 (0)1223 294048
E: techsupport@cytoell.com
W: www.oql.com

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
Symbols Glossary


EN ISO 15223-1:2021 - “Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements” (© International Organization for Standardization)		
Symbol	Title	Reference Number(s)
	en: Manufacturer	5.1.1
	en: Authorized representative in the European Community/European Union	5.1.2
	en: Use-by date	5.1.4
	en: Batch code	5.1.5
	en: Catalogue number	5.1.6
	en: Keep away from sunlight	5.3.2
	en: Temperature limit	5.3.7
	en: Consult instructions for use	5.4.3
	en: Consult electronic instructions for use	5.4.3
	en: Caution	5.4.4
	en: <i>In vitro</i> diagnostic medical device	5.5.1
	en: Contains sufficient for <n> tests	5.5.5

	en: Unique Device Identifier	5.7.10
EDMA symbols for IVD reagents and components, October 2009 revision		
Symbol	Title	Reference Number(s)
	en: Contents (or contains)	N/A

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Patents and Trademarks

Aquarius and Cytocell are registered trademarks of Cytocell Ltd.

Cytocell Ltd.
Oxford Gene Technology,
418 Cambridge Science Park,
Milton Road,
Cambridge, CB4 0PZ, UK
T: +44(0)1223 294048
E: probes@cytocell.com
W: www.cytocell.com