



A Sysmex Group Company

SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Canada Hazardous Products Act (HPA) and the Hazardous Products Regulation (HPR), as amended

Issuing Date 24-Oct-2025

Revision date 24-Oct-2025

Revision Number 1

1. Identification

Product identifier

Product Name Formamide

Other means of identification

Product Code(s) 770400-24/96, 770410-24/96, 770500-24/96, 770510-24/96 Universal Hybridization & Wash Kit.500075, 500082, 500083, 500084, 500085, 500086 SureSeq NGS Hybridization & Wash Kit.Complete Workflow Solutions (Universal)-78000*-24/96, 78010*-24/96, 780126-48, 780127-24/96, 79000*-24/96, 79010*-24/96, 890001-24/96

Synonyms None

Recommended use of the chemical and restrictions on use

Recommended use Additive
For research use only

Restrictions on use None known

Details of the supplier of the safety data sheet

Supplier Address

Oxford Gene Technology Inc.
(North America office)
520 White Plains Road, Suite
500
Tarrytown, NY 10591
USA
914 467 5285

Manufacturer Address

Oxford Gene Technology
Unit 5
4A Oxford Technology Park
Kidlington
OX5 1GN, United Kingdom
+44 (0)1865 856800
<http://www.ogt.com>

E-mail support@ogt.com

Emergency telephone number

Emergency telephone 914 467 5285

2. Hazard(s) identification

Classification of the substance or mixture

Reproductive toxicity	Category 1B
-----------------------	-------------

Label elements

**Danger****Hazard statements**

May damage fertility or the unborn child

Precautionary Statements - Prevention

Obtain special instructions before use

Do not handle until all safety precautions have been read and understood

Wear protective gloves, protective clothing, eye protection and face protection

Precautionary Statements - Response

IF exposed or concerned: Get medical advice/attention

Precautionary Statements - Storage

Store locked up

Precautionary Statements - Disposal

Dispose of contents and container in accordance with local, regional, national, and international regulations as applicable

Other information

No information available.

3. Composition/information on ingredients

Substance

Not applicable.

Mixture

Chemical name	CAS No.	Weight-%	Hazardous Material Information Review Act registry number (HMIRA registry #)	Date HMIRA filed and date exemption granted (if applicable)
Formamide	75-12-7	>99	-	

4. First-aid measures

Description of first aid measures**General advice**

Show this safety data sheet to the doctor in attendance.

Inhalation

Remove to fresh air.

Eye contact

Rinse thoroughly with plenty of water, also under the eyelids.

Skin contact

Wash with plenty of water.

Ingestion

Rinse mouth.

Most important symptoms and effects, both acute and delayed

Symptoms	None known.
Effects of Exposure	May cause adverse reproductive effects - such as birth defect, miscarriages, or infertility.

Indication of any immediate medical attention and special treatment needed

Note to physicians	Treat symptomatically.
---------------------------	------------------------

5. Fire-fighting measures

Suitable Extinguishing Media	Dry chemical, CO2, alcohol-resistant foam or water spray.
Unsuitable extinguishing media	None known.
Specific hazards arising from the chemical	No information available.
Hazardous combustion products	Toxic gases or vapors: Ammonia, Hydrogen cyanide, Carbon monoxide, Nitrogen oxides (NOx).
Explosion data	
Sensitivity to mechanical impact	None.
Sensitivity to static discharge	None.
Special protective equipment and precautions for fire-fighters	Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.

6. Accidental release measures**Personal precautions, protective equipment and emergency procedures**

Personal precautions	Ensure adequate ventilation.
-----------------------------	------------------------------

Methods and material for containment and cleaning up

Methods for containment	Prevent further leakage or spillage if safe to do so.
Methods for cleaning up	Dam up. Soak up with inert absorbent material. Pick up and transfer to properly labeled containers. Clean contaminated surface thoroughly.
Prevention of secondary hazards	Clean contaminated objects and areas thoroughly observing environmental regulations.

7. Handling and storage**Precautions for safe handling**

Advice on safe handling	Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes or clothing. Do not eat, drink or smoke when using this product. Remove contaminated clothing and shoes.
General hygiene considerations	Do not eat, drink or smoke when using this product. Wash hands before breaks and immediately after handling the product.

Conditions for safe storage, including any incompatibilities

Storage Conditions	Store at ambient conditions. Store locked up.
---------------------------	---

8. Exposure controls/personal protection

Control Parameters

Exposure Limits

Chemical name	Alberta	British Columbia	Ontario	Quebec
Formamide 75-12-7	TWA: 10 ppm; TWA: 18 mg/m ³ ; pSk	TWA: 10 ppm; Adverse reproductive effect Sk	TWA: 10 ppm; dSk	TWAEV: 1 ppm; Sd

Chemical name	Manitoba	New Brunswick	Newfoundland and Labrador	Nova Scotia
Formamide	TWA: 1 ppm; pSk	TWA: 10 ppm; pSk	TWA: 1 ppm; pSk	TWA: 1 ppm; pSk

Chemical name	Nunavut	Prince Edward Island	Saskatchewan	Yukon
Formamide	TWA: 10 ppm; STEL: 15 ppm; Sk	TWA: 1 ppm;	TWA: 10 ppm; STEL: 15 ppm; pSd	TWA: 20 ppm; TWA: 30 mg/m ³ ; STEL: 30 ppm; STEL: 45 mg/m ³ ;

Note See section 16 for terms and abbreviations.

Appropriate engineering controls

Engineering controls Showers
Eyewash stations
Ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/face protection Wear safety glasses with side shields (or goggles).

Hand protection Wear suitable gloves.

Skin and body protection Wear suitable protective clothing.

Respiratory protection Use appropriate respiratory protection. Consult with an industrial hygienist to determine the appropriate respiratory protection for your specific use of this material.

9. Physical and chemical properties

Information on basic physical and chemical properties

Appearance Colorless liquid
Physical state Liquid
Color Colorless
Odor Odorless
Odor threshold No information available

Property	Values	Remarks • Method
Melting point / freezing point	26 °C	
Initial boiling point and boiling range	218.3 °C	

Flammability		No data available
Flammability Limit in Air		
Upper flammability or explosive limits		No data available
Lower flammability or explosive limits		No data available
Flash point	152 °C	
Autoignition temperature	> 500 °C	
Decomposition temperature	> 140 °C	
SADT (°C)		No data available
pH		No data available
pH (as aqueous solution)		No data available
Kinematic viscosity		No data available
Dynamic viscosity	3.764 mPa s	@ 20 °C
Water solubility	Miscible in water	
Solubility(ies)		No data available
Partition Coefficient (n-octanol/water)	Log Kow: -0.82	@ 25 °C
Vapor pressure	0.06 hPa	@ 20 °C
Relative density		No data available
Bulk density		No data available
Liquid Density	1.13 g/cm ³ @ 20°C	No data available
Relative vapor density		No data available
Particle characteristics		No information available
Particle Size		No data available
Particle Size Distribution		No data available
Other information		
Molecular weight	No information available	
VOC content	No information available	
Softening point	No information available	

Information with regard to physical hazard classes

Explosives

Explosive properties Not an explosive.

Oxidizing properties

Not an oxidizer.

10. Stability and reactivity

Reactivity	None under normal use conditions.
Chemical stability	Stable under normal conditions.
Possibility of hazardous reactions	None under normal processing.
Conditions to avoid	Excessive heat. Extremes of temperature and direct sunlight.
Incompatible materials	Iodine, Pyridine, Sulfur trioxide, Natural rubber.
Hazardous decomposition products	None known based on information supplied.

11. Toxicological information

Information on likely routes of exposure

Product Information

Inhalation	Specific test data for the substance or mixture is not available.
Eye contact	Specific test data for the substance or mixture is not available.

Skin contact Specific test data for the substance or mixture is not available.

Ingestion Specific test data for the substance or mixture is not available.

Symptoms related to the physical, chemical and toxicological characteristics

Symptoms None known.

Acute toxicity No information available.

Numerical measures of toxicity

Component Information

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50
Formamide 75-12-7	= 5577 mg/kg (Rat)	= 6 g/kg (Rabbit)	> 21 mg/L (Rat) 4 h

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation No information available.

Serious eye damage/eye irritation No information available.

Respiratory or skin sensitization No information available.

Germ cell mutagenicity No information available.

Carcinogenicity Based on available data, the classification criteria are not met.

The table below indicates whether each agency has listed any ingredient as a carcinogen.

Chemical name	ACGIH	IARC	NTP	OSHA
Formamide 75-12-7	A3 - Confirmed animal carcinogen (with unknown relevance to humans)	-	-	-

Reproductive toxicity Contains a known or suspected reproductive toxin. Classification based on data available for ingredients. May damage fertility or the unborn child.

STOT - single exposure No information available.

STOT - repeated exposure No information available.

Aspiration hazard No information available.

12. Ecological information

Ecotoxicity Low toxicity to aquatic organisms. Based on available data, the classification criteria are not met.

Aquatic ecotoxicity

Component Information

Chemical name	Fish	Crustacea	Algae/aquatic plants	Toxicity to microorganisms
Formamide	LC50: =9135mg/L (96h,	EC50: >500mg/L (48h,	EC50: >500mg/L (72h,	-

	Brachydanio rerio)	Daphnia magna)	Desmodesmus subspicatus) EC50: >500mg/L (96h, Desmodesmus subspicatus)	
--	--------------------	----------------	--	--

Persistence and degradability Readily biodegradable.

Bioaccumulative potential

Chemical name	Partition coefficient	Bioconcentration factor (BCF)	Trophic magnification factor (TMF)
Formamide	-0.82	-	-

Mobility in soil No information available.

Other adverse effects No information available.

13. Disposal considerations

Disposal methods

Waste from residues/unused products Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.

Contaminated packaging Do not reuse empty containers.

14. Transport information

TDG Not regulated

IATA Not regulated

IMDG Not regulated

15. Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

International Regulations

The Montreal Protocol on Substances that Deplete the Ozone Layer Not applicable

The Stockholm Convention on Persistent Organic Pollutants Not applicable

The Rotterdam Convention Not applicable

International Inventories

Contact supplier for inventory compliance status

16. Other information

NFPA **Health hazards** 0 **Flammability** 1 **Instability** 0 **Special hazards** -
HMIS **Health hazards** * **Flammability** 1 **Physical hazards** 0 **Personal protection** X
Chronic Hazard Star Legend * = Chronic Health Hazard

Key or legend to abbreviations and acronyms used in the safety data sheet

List may include phrases which are not applicable to this product

ACGIH	American Conference of Governmental Industrial Hygienists
ADN	Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Europe)
ADR	Agreement concerning the International Carriage of Dangerous Goods by Road (Europe)
AIIC	Australian Inventory of Industrial Chemicals
ATE	Acute Toxicity Estimate
ASTM	American Society for the Testing of Materials
bar	Biological Reference Values for Chemical Compounds in the Work Area
BAT	Biological tolerance values for occupational exposure
BEL	Biological exposure limits
bw	Body weight
Ceiling	Maximum limit value
CMR	Carcinogen, Mutagen or Reproductive Toxicant
DOT	Department of Transportation (United States)
DSL	Domestic Substances List (Canada)
EmS	Emergency Schedule
ENCS	Existing and New Chemical Substances (Japan)
EPA	U.S. Environmental Protection Agency
GHS	Globally Harmonized System
HMIS	Hazardous Materials Identification System
IARC	International Agency for Research on Cancer
IATA	International Air Transport Association
IBC	International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk
ICAO	International Civil Aviation Organization
IECSC	Inventory of Existing Chemical Substances in China
IMDG	International Maritime Dangerous Goods
IMO	International Maritime Organization
ISO	International Organization for Standardization
KECI	Korean Existing Chemicals Inventory
LC50	Lethal Concentration to 50% of a test population
LD50	Lethal Dose to 50% of a test population (Median Lethal Dose)
MARPOL	International Convention for the Prevention of Pollution from Ships
NFPA	National Fire Protection Association
n.o.s.	Not Otherwise Specified
NOAEC	No Observed Adverse Effect Concentration
NOAEL	No Observed Adverse Effect Level
NOELR	No Observable Effect Loading Rate
NZIoC	New Zealand Inventory of Chemicals
OECD	Organization for Economic Cooperation and Development
OEL	Occupational exposure limits
PBT	Persistent, Bioaccumulative and Toxic substance
PICCS	Philippines Inventory of Chemicals and Chemical Substances
PMT	Persistent, Mobile and Toxic
PPE	Personal protective equipment
QSAR	Quantitative Structure Activity Relationship
RID	Agreement concerning the International Carriage of Dangerous Goods by Rail (Europe)
SADT	Self-Accelerating Decomposition Temperature
SAR	Structure-activity relationship

SDS	Safety Data Sheet
SL	Surface Limit
STEL	Short Term Exposure Limit
STOT RE	Specific target organ toxicity - Repeated exposure
STOT SE	Specific target organ toxicity - Single exposure
TCSI	Taiwan Chemical Substance Inventory
TDG	Transport of Dangerous Goods (Canada)
TSCA	Toxic Substances Control Act (United States)
TWA	Time-Weighted Average
UN	United Nations
VOC	Volatile organic compounds
vPvB	Very Persistent and Very Bioaccumulative
vPvM	Very Persistent and Very Mobile
As	Allergenic substance
DS	Dermal Sensitizer
Ot	Ototoxicant
pOt	Ototoxicant - potential to cause hearing disorders
PS	Photosensitizer
RS	Respiratory Sensitizer
S	Sensitizer
poS	Sensitizer - capable of causing occupational asthma
Sa	Simple asphyxiant
Sd	Skin designation
pSd	Skin designation - potential for cutaneous absorption
Sdv	Skin designation - vacated
Sk	Skin notation
dSk	Skin notation - danger of cutaneous absorption
pSk	Skin notation - potential for cutaneous absorption

Key literature references and sources for data used to compile the SDS

U.S. Environmental Protection Agency ChemView Database
 European Food Safety Authority (EFSA)
 U.S. Environmental Protection Agency
 Acute Exposure Guideline Level(s) (AEGl(s))
 U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act
 U.S. Environmental Protection Agency High Production Volume Chemicals
 Food Research Journal
 Hazardous Substance Database
 International Uniform Chemical Information Database (IUCLID)
 Japan GHS Classification
 Australia National Industrial Chemicals Notification and Assessment Scheme (NICNAS)
 NIOSH (National Institute for Occupational Safety and Health)
 National Library of Medicine's ChemID Plus (NLM CIP)
 National Library of Medicine's PubMed database (NLM PUBMED)
 U.S. National Toxicology Program (NTP)
 New Zealand's Chemical Classification and Information Database (CCID)
 International Organization for Economic Co-operation and Development (OECD) Environment, Health, and Safety Publications
 International Organization for Economic Co-operation and Development (OECD) High Production Volume Chemicals Program
 International Organization for Economic Co-operation and Development (OECD) Screening Information Data Set
 United Nations World Health Organization (WHO)

Issuing Date 24-Oct-2025

Revision date 24-Oct-2025

Revision Note Initial Release.

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

End of Safety Data Sheet