

Quality Control Test Certificate

Product name	SureSeq FFPE DNA Repair Mix
Product no.	500079
Batch no.	10011027
Certificate issue no.	1

Acceptance Criteria

The new stock of SureSeq FFPE DNA Repair Mix was tested using a SureSeq NGS Cancer Panel for target enrichment with a human genomic DNA library. Treated and untreated DNA was run on an Illumina sequencer to verify that there was an improvement after treatment and met the below series of quality specifications.

FFPE DNA Repair Mix	Specification	Pass/Fail
Post-PCR1 Concentration Improvement:	>10%	PASS
% On Target Bases Improvement:	>10%	PASS
Mean Target Coverage Improvement:	>10%	PASS

This product has been approved for release by:

Name: Vincent Wong Position: Quality Associate Date: 17 August 2018

This product has been manufactured and tested by Oxford Gene Technology (Operations) Ltd. under an ISO 9001:2015 certified management system.

Oxford Gene Technology

T: +44(0)1865 856826 (US: 914-467-5285) E: products@ogt.com W: www.ogt.com
 Technical Support: support@ogt.com

The purchaser has the non-transferable right to use and consume the product for RESEARCH USE ONLY AND NOT for DIAGNOSTICS PROCEDURES. It is not intended for use, and should not be used for the diagnosis, prevention, monitoring, treatment or alleviation of any disease or condition, or for the investigation of any physiological process, in any identifiable human, or for any other medical purpose.