

## QC Certificate

Product Name	SureSeq NGS Hyb & Wash Kit
Product Code	500075
Lot Number	079907

### Acceptance Criteria

*The SureSeq NGS Hyb & Wash Kit was used with the short and long hyb. protocols alongside a SureSeq Core MPN and Custom Cancer MyPanel respectively; and then sequenced on an Illumina sequencer to verify that it meets the below series of quality specifications:*

**Number of Lanes:** 16  
**Starting DNA Input:** 1000ng

		Specification	Pass/Fail
Mean Target Coverage:	Core MPN:	> 3000x	Pass
	Custom MyPanel:	> 1000x	Pass
% Target Bases Covered at >20% Mean:		>99.0%	Pass

Name	Position	Date
Kaajal Reeves	QC Manager	29 April 2022

This product has been released by Oxford Gene Technology (Operations) Ltd. under an ISO 9001:2015 certified management system.

### OGT

T: +44(0)1865 856800 (US: 914-467-5285) E: [products@ogt.com](mailto:products@ogt.com) W: [www.ogt.com](http://www.ogt.com)

Technical Support: [support@ogt.com](mailto:support@ogt.com)

The purchaser has the non-transferable right to use and consume the product for RESEARCH USE ONLY AND NOT for DIAGNOSTICS PROCEDURES. It is not intended for use, and should not be used for the diagnosis, prevention, monitoring, treatment or alleviation of any disease or condition, or for the investigation of any physiological process, in any identifiable human, or for any other medical purpose.