

## QC Certificate

|              |                        |
|--------------|------------------------|
| Product Name | Breast Cancer Plus CNV |
| Product Code | 602028                 |
| Lot Number   | 486280741              |

### Acceptance Criteria

**Number of Lanes:** 16  
**Starting Input DNA:** 1000 ng

| Breast Cancer + CNV Panel       | Specification | Pass/Fail |
|---------------------------------|---------------|-----------|
| Mean Target Coverage:           | ≥1000x        | Pass      |
| ≥98.0% Target Bases Covered at: | ≥500x         | Pass      |
| % Reads Covered at >20% Mean:   | ≥98.5%        | Pass      |

| Name          | Position   | Date        |
|---------------|------------|-------------|
| Kaajal Reeves | QC Manager | 27 May 2022 |

This product has been released by Oxford Gene Technology (Operations) Ltd. under an ISO 9001:2015 certified management system.

### Oxford Gene Technology

T: +44(0)1865 856800 (US: 914-467-5285) E: [products@ogt.com](mailto:products@ogt.com) W: [www.ogt.com](http://www.ogt.com)

Technical Support: [support@ogt.com](mailto:support@ogt.com)

The purchaser has the non-transferable right to use and consume the product for RESEARCH USE ONLY AND NOT for DIAGNOSTICS PROCEDURES. It is not intended for use, and should not be used for the diagnosis, prevention, monitoring, treatment or alleviation of any disease or condition, or for the investigation of any physiological process, in any identifiable human, or for any other medical purpose.