



SAFETY DATA SHEET

This safety data sheet was created pursuant to the requirements of:
Canada Hazardous Products Act (HPA) and the Hazardous Products Regulation (HPR), as
amended

Issuing Date 22-Jan-2024

Revision date 23-Mar-2026

Revision Number 2

1. Identification

Product identifier

Product Name HEAT PRETREATMENT SOLUTION PH 7

Other means of identification

Product Code(s) LPS100A

Synonyms None

Recommended use of the chemical and restrictions on use

Recommended use Laboratory chemicals
For professional use only

Restrictions on use None known

Details of the supplier of the safety data sheet

Supplier Address

Oxford Gene Technology Inc.
(North America office)
520 White Plains Road, Suite
500
Tarrytown, NY 10591
USA
914 467 5285

Manufacturer Address

CytoCell Ltd., Oxford Gene Technology
418 Cambridge Science Park, Milton Road,
Cambridge
CB4 0PZ, United Kingdom
T: +44 (0)1223 294048
F: +44 (0)1223 294986
probes@cytoCell.com
<http://www.ogt.com>

E-mail support@ogt.com

Emergency telephone number

Emergency telephone 914 467 5285

2. Hazard(s) identification

Classification of the substance or mixture

| | |
|--|-------------|
| Serious eye damage/eye irritation | Category 2A |
| Specific target organ toxicity (repeated exposure) | Category 2 |

Label elements

**Warning****Hazard statements**

Causes serious eye irritation
May cause damage to organs through prolonged or repeated exposure

Precautionary Statements - Prevention

Wash hands thoroughly after handling
Do not breathe dust, fume, gas, mist, vapors and spray
Wear protective gloves, protective clothing, eye protection and face protection

Precautionary Statements - Response

Get medical advice/attention if you feel unwell

Eyes

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
If eye irritation persists: Get medical advice and attention

Precautionary Statements - Disposal

Dispose of contents and container in accordance with local, regional, national, and international regulations as applicable

Other information

The product does not contain any substance(s) classified as PBT or vPvB.

3. Composition/information on ingredients

Substance

Not applicable.

Mixture

| Chemical name | CAS No. | Weight-% | Hazardous Material Information Review Act registry number (HMIRA registry #) | Date HMIRA filed and date exemption granted (if applicable) |
|--|-------------|----------|--|---|
| Ethylenediaminetetraacetic acid tetrasodium salt hydrate | 194491-31-1 | 1 - 3 | - | |

4. First-aid measures

Description of first aid measures**General advice**

Show this safety data sheet to the doctor in attendance.

Inhalation

Remove person to fresh air and keep comfortable for breathing. Remove to fresh air.

Eye contact

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Keep eye wide open while rinsing. Do not rub affected area. Get medical attention if irritation develops and persists.

| | |
|---|---|
| Skin contact | Wash with plenty of water. |
| Ingestion | Rinse mouth. Never give anything by mouth to an unconscious person. Do NOT induce vomiting. Call a physician. |
| Self-protection of the first aider | Avoid contact with skin, eyes or clothing. Wear personal protective clothing (see section 8). |

Most important symptoms and effects, both acute and delayed

| | |
|----------------------------|--|
| Symptoms | May cause redness and tearing of the eyes. Burning sensation. |
| Effects of Exposure | May cause damage to organs through prolonged or repeated exposure. |

Indication of any immediate medical attention and special treatment needed

| | |
|---------------------------|------------------------|
| Note to physicians | Treat symptomatically. |
|---------------------------|------------------------|

5. Fire-fighting measures

| | |
|---|--|
| Suitable Extinguishing Media | Dry chemical, CO2, water spray or regular foam. |
| Unsuitable extinguishing media | Do not use a solid water stream as it may scatter and spread fire. |
| Specific hazards arising from the chemical | None known based on information supplied. |
| Hazardous combustion products | Carbon oxides. |
| Explosion data | |
| Sensitivity to mechanical impact | None. |
| Sensitivity to static discharge | None. |
| Special protective equipment and precautions for fire-fighters | Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protective equipment. |

6. Accidental release measures**Personal precautions, protective equipment and emergency procedures**

| | |
|-----------------------------|---|
| Personal precautions | Avoid contact with skin, eyes or clothing. Use personal protective equipment as required. |
| Other information | Refer to protective measures listed in Sections 7 and 8. |

Methods and material for containment and cleaning up

| | |
|--|--|
| Methods for containment | Prevent further leakage or spillage if safe to do so. |
| Methods for cleaning up | Dam up. Soak up with inert absorbent material. Pick up and transfer to properly labeled containers. Clean contaminated surface thoroughly. |
| Prevention of secondary hazards | Clean contaminated objects and areas thoroughly observing environmental regulations. |

7. Handling and storage**Precautions for safe handling**

| | |
|--------------------------------|--|
| Advice on safe handling | Wear personal protective clothing (see section 8). Wash hands thoroughly after handling. Handle in accordance with good industrial hygiene and safety practice. Avoid contact with |
|--------------------------------|--|

skin, eyes or clothing. Do not eat, drink or smoke when using this product.

General hygiene considerations Avoid contact with skin, eyes or clothing. Wear suitable gloves and eye/face protection. Do not eat, drink or smoke when using this product.

Conditions for safe storage, including any incompatibilities

Storage Conditions Keep containers tightly closed in a dry, cool and well-ventilated place.

8. Exposure controls/personal protection

Control Parameters

Exposure Limits This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

Appropriate engineering controls

Engineering controls Showers
Eyewash stations
Ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/face protection Handling of larger amounts: Wear protective eye glasses for protection against liquid splashes. Wear safety glasses with side shields (or goggles).

Hand protection Wear protective gloves. To protect the wearer, gloves must be the correct fit and be used properly. Ensure that the breakthrough time of the glove material is not exceeded. Refer to glove supplier for information on breakthrough time for specific gloves. Wear suitable gloves.

Skin and body protection Wear suitable protective clothing.

Respiratory protection Use appropriate respiratory protection. Consult with an industrial hygienist to determine the appropriate respiratory protection for your specific use of this material.

9. Physical and chemical properties

Information on basic physical and chemical properties

Appearance Colorless liquid
Physical state Liquid
Color Colorless
Odor Characteristic
Odor threshold No information available

| <u>Property</u> | <u>Values</u> | <u>Remarks • Method</u> |
|--|---------------|-------------------------|
| Melting point / freezing point | | No data available |
| Initial boiling point and boiling range | | No data available |
| Flammability | | No data available |
| Flammability Limit in Air | | |
| Upper flammability or explosive limits | | Not applicable |
| Lower flammability or explosive limits | | Not applicable |
| Flash point | | No data available |
| Autoignition temperature | | No data available |

| | | |
|--|-------|--------------------------|
| Decomposition temperature | | No data available |
| SADT (°C) | | No data available |
| pH | 6 - 8 | |
| pH (as aqueous solution) | | No data available |
| Kinematic viscosity | | No data available |
| Dynamic viscosity | | No data available |
| Water solubility | | No data available |
| Solubility(ies) | | No data available |
| Partition Coefficient (n-octanol/water) | | No data available |
| Vapor pressure | | No data available |
| Relative density | | No data available |
| Bulk density | | No data available |
| Liquid Density | | No data available |
| Relative vapor density | | No data available |
| Particle characteristics | | No information available |
| Particle Size | | No data available |
| Particle Size Distribution | | No data available |

Other information

| | |
|------------------|---------------------------|
| Molecular weight | No information available |
| VOC content | No information available. |
| Softening point | No information available |

Information with regard to physical hazard classes**Explosives**

Explosive properties No information available.

Oxidizing properties

No information available.

10. Stability and reactivity

| | |
|------------------------------------|---|
| Reactivity | None under normal use conditions. |
| Chemical stability | Stable under normal conditions. |
| Possibility of hazardous reactions | None under normal processing. |
| Conditions to avoid | Protect from direct sunlight. Extremely high or low temperatures. |
| Incompatible materials | None known based on information supplied. |
| Hazardous decomposition products | Carbon oxides. |

11. Toxicological information**Information on likely routes of exposure****Product Information**

| | |
|---------------------|---|
| Inhalation | Specific test data for the substance or mixture is not available. May cause irritation of respiratory tract. |
| Eye contact | Specific test data for the substance or mixture is not available. Causes serious eye irritation. (based on components). May cause redness, itching, and pain. |
| Skin contact | Specific test data for the substance or mixture is not available. May cause irritation. Prolonged contact may cause redness and irritation. |
| Ingestion | Specific test data for the substance or mixture is not available. Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea. |

Symptoms related to the physical, chemical and toxicological characteristics

| | |
|---------------------------------------|---|
| Symptoms | May cause redness and tearing of the eyes. Burning sensation. |
| <u>Acute toxicity</u> | No information available. |
| Numerical measures of toxicity | No information available. |

Delayed and immediate effects as well as chronic effects from short and long-term exposure

| | |
|--|--|
| Skin corrosion/irritation | No information available. |
| Serious eye damage/eye irritation | Classification based on data available for ingredients. Causes serious eye irritation. |
| Respiratory or skin sensitization | No information available. |
| Germ cell mutagenicity | No information available. |
| Carcinogenicity | No information available. |
| Reproductive toxicity | No information available. |
| STOT - single exposure | No information available. |
| STOT - repeated exposure | May cause damage to organs through prolonged or repeated exposure. |
| Aspiration hazard | No information available. |

12. Ecological information

Ecotoxicity The environmental impact of this product has not been fully investigated.

Persistence and degradability No information available.

Bioaccumulative potential No information available.

Mobility in soil No information available.

Other adverse effects No information available.

13. Disposal considerations**Disposal methods**

Waste from residues/unused products Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.

Contaminated packaging Do not reuse empty containers.

14. Transport information

TDG Not regulated

IATA Not regulated

IMDG Not regulated

15. Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

International Regulations

The Montreal Protocol on Substances that Deplete the Ozone Layer Not applicable

The Stockholm Convention on Persistent Organic Pollutants Not applicable

The Rotterdam Convention Not applicable

International Inventories

Contact supplier for inventory compliance status

16. Other information

| | | | | |
|-----------------------------------|---------------------------|-----------------------|---------------------------|------------------------------|
| NFPA | Health hazards 2 | Flammability 0 | Instability 0 | Special hazards - |
| HMIS | Health hazards 2 * | Flammability 0 | Physical hazards 0 | Personal protection X |
| Chronic Hazard Star Legend | * = Chronic Health Hazard | | | |

Key or legend to abbreviations and acronyms used in the safety data sheet

List may include phrases which are not applicable to this product

| | |
|---------|---|
| ACGIH | American Conference of Governmental Industrial Hygienists |
| ADN | Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Europe) |
| ADR | Agreement concerning the International Carriage of Dangerous Goods by Road (Europe) |
| AIIC | Australian Inventory of Industrial Chemicals |
| ATE | Acute Toxicity Estimate |
| ASTM | American Society for the Testing of Materials |
| bar | Biological Reference Values for Chemical Compounds in the Work Area |
| BAT | Biological tolerance values for occupational exposure |
| BEL | Biological exposure limits |
| bw | Body weight |
| Ceiling | Maximum limit value |
| CMR | Carcinogen, Mutagen or Reproductive Toxicant |
| DOT | Department of Transportation (United States) |
| DSL | Domestic Substances List (Canada) |
| EmS | Emergency Schedule |
| ENCS | Existing and New Chemical Substances (Japan) |
| EPA | U.S. Environmental Protection Agency |
| GHS | Globally Harmonized System |
| HMIS | Hazardous Materials Identification System |
| IARC | International Agency for Research on Cancer |
| IATA | International Air Transport Association |
| IBC | International Code for the Construction and Equipment of Ships carrying Dangerous |

| | |
|---------|---|
| | Chemicals in Bulk |
| ICAO | International Civil Aviation Organization |
| IECSC | Inventory of Existing Chemical Substances in China |
| IMDG | International Maritime Dangerous Goods |
| IMO | International Maritime Organization |
| ISO | International Organization for Standardization |
| KECI | Korean Existing Chemicals Inventory |
| LC50 | Lethal Concentration to 50% of a test population |
| LD50 | Lethal Dose to 50% of a test population (Median Lethal Dose) |
| MARPOL | International Convention for the Prevention of Pollution from Ships |
| NFPA | National Fire Protection Association |
| n.o.s. | Not Otherwise Specified |
| NOAEC | No Observed Adverse Effect Concentration |
| NOAEL | No Observed Adverse Effect Level |
| NOELR | No Observable Effect Loading Rate |
| NZIoC | New Zealand Inventory of Chemicals |
| OECD | Organization for Economic Cooperation and Development |
| OEL | Occupational exposure limits |
| PBT | Persistent, Bioaccumulative and Toxic substance |
| PICCS | Philippines Inventory of Chemicals and Chemical Substances |
| PMT | Persistent, Mobile and Toxic |
| PPE | Personal protective equipment |
| QSAR | Quantitative Structure Activity Relationship |
| RID | Agreement concerning the International Carriage of Dangerous Goods by Rail (Europe) |
| SADT | Self-Accelerating Decomposition Temperature |
| SAR | Structure-activity relationship |
| SDS | Safety Data Sheet |
| SL | Surface Limit |
| STEL | Short Term Exposure Limit |
| STOT RE | Specific target organ toxicity - Repeated exposure |
| STOT SE | Specific target organ toxicity - Single exposure |
| TCSI | Taiwan Chemical Substance Inventory |
| TDG | Transport of Dangerous Goods (Canada) |
| TSCA | Toxic Substances Control Act (United States) |
| TWA | Time-Weighted Average |
| UN | United Nations |
| VOC | Volatile organic compounds |
| vPvB | Very Persistent and Very Bioaccumulative |
| vPvM | Very Persistent and Very Mobile |
| As | Allergenic substance |
| DS | Dermal Sensitizer |
| Ot | Ototoxicant |
| pOt | Ototoxicant - potential to cause hearing disorders |
| PS | Photosensitizer |
| RS | Respiratory Sensitizer |
| S | Sensitizer |
| poS | Sensitizer - capable of causing occupational asthma |
| Sa | Simple asphyxiant |
| Sd | Skin designation |
| pSd | Skin designation - potential for cutaneous absorption |
| Sdv | Skin designation - vacated |
| Sk | Skin notation |
| dSk | Skin notation - danger of cutaneous absorption |
| pSk | Skin notation - potential for cutaneous absorption |

Key literature references and sources for data used to compile the SDS

U.S. Environmental Protection Agency ChemView Database
European Food Safety Authority (EFSA)
U.S. Environmental Protection Agency
Acute Exposure Guideline Level(s) (AEGL(s))
U.S. Environmental Protection Agency Federal Insecticide, Fungicide, and Rodenticide Act
U.S. Environmental Protection Agency High Production Volume Chemicals
Food Research Journal
Hazardous Substance Database
International Uniform Chemical Information Database (IUCLID)
Japan GHS Classification
Australia National Industrial Chemicals Notification and Assessment Scheme (NICNAS)
NIOSH (National Institute for Occupational Safety and Health)
National Library of Medicine's ChemID Plus (NLM CIP)
National Library of Medicine's PubMed database (NLM PUBMED)
U.S. National Toxicology Program (NTP)
New Zealand's Chemical Classification and Information Database (CCID)
International Organization for Economic Co-operation and Development (OECD) Environment, Health, and Safety Publications
International Organization for Economic Co-operation and Development (OECD) High Production Volume Chemicals Program
International Organization for Economic Co-operation and Development (OECD) Screening Information Data Set
United Nations World Health Organization (WHO)

Issuing Date 22-Jan-2024

Revision date 23-Mar-2026

Revision Note Change to composition. Change in the mixture classification.

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

End of Safety Data Sheet